

Product datasheet for **SC332784**

CD229 (LY9) (NM_001261457) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: CD229 (LY9) (NM_001261457) Human Untagged Clone
Tag: Tag Free
Symbol: CD229
Synonyms: CD229; hly9; mLY9; SLAMF3
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC332784 representing NM_001261457.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGTGGCACAAAGAGTCACACAGATGACTGGGCTCCTGGCCTTTCTCCAGTAAGCCACAGAGGAGT
CAGCTGCAAATATTCTCTTCTGTCTACAGACCTCTCTCCTTCTCCTGCTCATGGGACTAAGAGCCTCT
GGAAAGGACTCAGCCCCAACAGTGGTGTGAGGGATCCTAGGGGGTTCGTGACTCTCCCCCTAACATC
TCAGTAGACACAGAGATTGAGAACGTCTCTGGATTGGTCCAAAAATGCTCTTGCTTTCCGACGTCCC
AAAGAAAATGTAACCATTATGGTCAAAGCTACCTGGGCCGACTAGACATACCAAGTGGAGTTACTCC
CTGTGCATCAGCAATCTGACTCTGAATGATGCAGGATCCTACAAAGCCCAGATAAAACCAAGGAATTTT
GAAGTCACCACTGAGGAGGAATTCACCCTGTTCGTCTATGAGCAGCTGCAGGAGCCCAAGTCACCATG
AAGTCTGTGAAGGTGCTGAGAATTCTCTGTAACATCACTCTAATGTGCTCCGTGAAGGGGGCAGAG
AAAAGTGTCTGTACAGCTGGACCCCAAGGGAACCCCATGCTTCTGAGTCCAATGGAGGCTCCATTTCT
ACCGTCTCCCGAACCATGTGACCCAGACCTGCCATACATCTGCACAGCCAGAACCCCGTCAGCCAG
AGAAGTCCCTCCCTGTCCATGTTGGGCAGTCTGTACAGATCCAGGAGCCTCCAGAGGAGGAACAACG
GGGAGACTGTGGTAGGGGTCTGGGAGAGCCAGTCACCCTGCCACTTGCCTCCAGCCTGCCGGGAC
ACAGAGAAGGTTGTCTGGTTGTTAACACATCCATCATTAGCAAAGAGAGGGAAGAAGCAGCAACGGCA
GATCCACTCATTAAATCCAGGGATCCTTACAAGAACAGGGTGTGGGTCTCCAGCCAGGACTGCTCCCTG
AAGATCAGCCAGCTGAAGATAGAGGACGCGGCCCTACCATGCCTACGTGTGCTCAGAGGCTCCAGC
GTCACCAGCATGACACATGTCACCCTGCTCATCTACCGACTGAGAGAAAACAAAGCTTTGGATTGGG
TTGTTCTGATGGTTTGCCTTCTGTGCGTTGGGATCTTCAGTGGTGCATTTGGAAGCGAAAAGGACGG
TGTTTCAGTCCCAGCCTTCTGTTCCAGCCAAGCTGAGGCCCCAGCGGATACACCAGAACCCACAGCTGGC
CACACGCTATACTCTGTGCTCTCCCAAGGATATGAGAAGCTGGACTCCCTCAGGCCTGCCAGGCAA
CAGCCTACACCCACCTCAGACAGCAGCTCTGACAGCAACCTCACAAGTGGAGGATGAGGACAGGCCT
GAGGTGCACAAGCCCATCAGTGAAGATATGAGGTATTTGACCAGGTCCTCAGGAGGGCGCTGGACAT
GACCCAGCCCCTGAGGGCAAGCAGACTATGATCCCGTCACTCCATATGTCACGGAAGTTGAGTCTGTG
GTTGGAGAGAACCACATGTATGCACAAGTGTCAACTACAGGGAAGACCCAGTTTCTCAGAAGGAA
GAGAGCTCAGCCACAATCTACTGCTCCATACGGAACCTCAGGTGGTGCCACCACCACAACAAGATGAT
CTTGAGATTCTGAAAGTCTACCTATGAAAATTCACCTGA
  
```

Restriction Sites: SgfI-MluI
ACCN: NM_001261457



[View online »](#)

Insert Size:	1698 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001261457.1</u>
RefSeq Size:	2280 bp
RefSeq ORF:	1698 bp
Locus ID:	4063
Cytogenetics:	1q23.3
Protein Families:	Transmembrane
MW:	62.1 kDa
Gene Summary:	<p>LY9 belongs to the SLAM family of immunomodulatory receptors (see SLAMF1; MIM 603492) and interacts with the adaptor molecule SAP (SH2D1A; MIM 300490) (Graham et al., 2006 [PubMed 16365421]).[supplied by OMIM, Mar 2008]</p> <p>Transcript Variant: This variant (4) lacks an exon in the coding region, but maintains the reading frame, compared to variant 1. The encoded isoform (d) is shorter than isoform a.</p>