

## Product datasheet for **SC332805**

### **KIF18B (NM\_001265577) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	KIF18B (NM_001265577) Human Untagged Clone
Tag:	Tag Free
Symbol:	KIF18B
Vector:	pCMV6-Entry (PS100001)



[View online »](#)

**Fully Sequenced ORF:** >SC332805 representing NM\_001265577.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGCAGTGGAGGACAGCACGCTGCAAGTAGTGGTACGGGTGCGGCCCCACCCTCGGGAGCTGGAC
AGTCAGCGGGCCAGTGGTTACAGTGGTGGACGAGCGGTGCTGGTGTTTAACCTGAGGAGCCCGAT
GGAGGGTTCCTGGCCTGAAATGGGGTGGCACCATGATGGCCCAAGAAGAAGGGCAAAGACCTGACG
TTTGTCTTTGACCGGGTCTTTGGCGAGGGGCCACCCAACAGGACGTGTTCCAGCACACCACGCACAGC
GTCCTGGACAGCTTCTCCAGGGCTACAACCTGCTCAGTGTTCCTACGGGGCCACCGGGGCTGGGAAG
ACACACACCATGCTGGGAAGGGAGGGGGACCCCGCATCATGTACCTGACCACCTGGAAGTGTACAGG
CGCCTGGAGGCCCGCAGCAGGAGAAGCACTTCGAGGTGCTCATCAGCTACCAGGAGGTGTATAATGAA
CAGATCCATGACCTCCTGGAGCCCAAGGGGCCCTTGCCATCCGCGAGGACCCCGACAAGGGGGTGGTG
GTGCAAGGACTTTCTTTCCACCAGCCAGCCTCAGCCGAGCAGCTGCTGGAGATACTGACCAGGGGGAAC
CGTAACCCGACGCAGCACCCACTGATGCCAACCGGACTTCTCCCGCTCCCATGCCATCTTCAGATC
TTTGTGAAGCAGCAGGACCGGGTCCAGGACTGACCCAGGCTGTCCAGGTGGCCAAGATGAGCCTGATT
GACCTGGCTGGCTCAGAGCGGGCATCCAGCACCCATGCGAAGGGGGAGCGGCTGCGGGAGGGGGCCAAC
ATCAACCCTCTCTGCTGGCGCTCATCAACGTCTCAATGCCTTGCCGATGCAAAGGGCCGCAAGACC
CATGTGCCCTACCGGGACAGCAAAGTACCCGCTGCTCAAAGACTCCCTCGGGGGCAACTGCCGACA
GTGATGATCGCTGCCATCAGCCCTCCAGCCTGACCTACGAGGACACGTACAACACCCTCAAATATGCC
GACCGGGCCAAGGAGATCAGGCTCTCGTGAAAGAGCAATGTGACCAGCCTGGACTGTCACATCAGCCAG
TATGCTACCATCTGCCAACAGCTCCAGGCTGAGGTAGCCGCTCTGAGGAAGAAGTCCAAGTGTATGAG
GGGGGAGGCCAGCCCCACCACAGGACCTCCAGGATCTCCAAGTCGGGACCACCACAGAACACCAG
CCCTGCACCCAGAGCTCCCTGCAGGGCCTAGAGCCCTTCAAGAGGAGAGTCTGGGGATGGAGGCCCA
GTGGAGAGGGCCATGGAAGGGAAGTCTTACAGACCAGGAGCAGTCCCAGAGGATGAGGATGAAGGCCCA
GCTGAGGAGGTTCCAACCCAGATGCCAGAGCAGAACCCACACATGCACTGCCAGATCCCTCGCCTG
ACCCTGCAGCCCAAGCCAGTCTGTTGGGCCACTTCTCAGCACGGGAAGTGGATGGGGACCGTTCTAAGCAG
TTGGCCCTAAAGGTGCTGTGCGTTGCCAGCGGAGTACTCCCTGCTCCAAGCAGCCAACCTCCTGACG
CCCACATGATCACAGAGTTTGGAGCCCTACAGCAGCTGGTGCAAGAGGAAAAAATTGAGCCTGGGGCA
GAGGCCTTGAGGACTTCAAGCCTGGCCAGGGGGCACCTCGGCTCAGGAGCTGTGTTTCAAGTCAATC
CCTGTGCCGTCTCTCTGCCCCAGAGCCTCCAGGATACACTGGCCCTGTGACCCGGACTATGGCGAGG
CGACTGAGTGGCCCCCTGCACACCCTGGGAATCCCGCTGGACCCAAGTGCACCCAGCCAGGGGTCC
CGATGGCCCATGGAGAAGAAGAGGAGGAGACCAAGCGCCTTGGAGGCAGACAGTCCCATGGCCCCAAG
CGGGGCACCAAGCGCCAGCGCCAGTCTTCTGCCCTGCCTAAGGAGAGGGTCTCTGCTGACACCCAA
CCTTACAGGGGGCCAGCACCCCAAGGAGAAAGGGCCTCCTCCCCCTGCCATTCCCCTCGCGTTTGC
CCAGCCACAGTCAAAAAGCCGGGTGCCCTGGGCCCTTCCGCCATGCAGAAGTGTCCACCCCGCTG
GCTCTGCCACTCGAGACCTCAATGCCACCTTTGATCTCTCTGAGGAGCCTCCCTCAAAGCCAGTTTC
CATGAATGCATTGGCTGGGACAAAATACCCAGGAGCTGAGCAGGCTGGACCAGCCCTTCAATCCCCAGG
GCACCTGTGCCCTGTTACCATGAAGGGCCCCAAGCCAACATCTTCCCTCCCTGGGACCTTGCCTGC
AAGAAGAAGCGCGTTGCGAGTTCTCAGTCTCCATGGCCGACGCGCATCGCCCGCTCCCCAGCAGC
ACTTTGAAGAGGCCAGCTGGGCCCTTGTACTCCAGAGCTGCCCTTGTGTCCTGTGCCCTAGCAAC
CGGAGGAATGAAAGGACCTCATCAGGGTGGGAGAGCACTCTCAGCAGGGAACGGCGTACCAAGGTG
TCCTGA
  
```

**Restriction Sites:** SgfI-MluI

**ACCN:** NM\_001265577

**Insert Size:** 2559 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001265577.1</a></u>
<b>RefSeq Size:</b>	4045 bp
<b>RefSeq ORF:</b>	2559 bp
<b>Locus ID:</b>	146909
<b>UniProt ID:</b>	<u><a href="#">Q86Y91</a></u>
<b>Cytogenetics:</b>	17q21.31
<b>Protein Families:</b>	Druggable Genome
<b>MW:</b>	93 kDa
<b>Gene Summary:</b>	<p>In complex with KIF2C, constitutes the major microtubule plus-end depolymerizing activity in mitotic cells. Its major role may be to transport KIF2C and/or MAPRE1 along microtubules. [UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) encodes the longer isoform (1).</p>