

Product datasheet for **SC332811**

RTN3 (NM_001265590) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RTN3 (NM_001265590) Human Untagged Clone
Tag:	Tag Free
Symbol:	RTN3
Synonyms:	ASYIP; HAP; NSPL2; NSPLII; RTN3-A1
Vector:	pCMV6-Entry (PS100001)



[View online »](#)

Fully Sequenced ORF: >SC332811 representing NM_001265590.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGCGGAGCCGTCGGCGGCCACTCAGTCCCATTCCATCTCTCGTCGTCTTCGGAGCCGAGCCGTCC
GCGCCCGCGCGCGGGAGCCAGGAGCCTGCCCGCCCTGGGGACGAAGAGCTGCAGCTCCTCCTGT
GCGGCCAGTTTCCAGAGCATCCTGCTTTTCTCTCAAAGAAAATGGTCAAGTGGAAAGCAAATAGAT
AAAGAGACCAAGAACCCAAATGGGGTATCAAGTAGGGAGGCTAAAAGTGCATTGGATGCTGATGACAGA
TTCACTTTGCTGACAGCCAGAAAACCCTACTGAGTACTCTAAGGTAGAAGGCATTTATACATATTCT
TTGTCTCCATCCAAAGTTTCAGGAGATGATGTTATTGAAAAGGATTCCCTGAATCACCATTGGAAGTA
ATTATTGACAAAGCAGCATTTGACAAAGAATTTAAAGACTCATATAAGGAGAGCACAGATGATTTTGGT
AGCTGGTCTGTGCACACTGATAAAGAATCATCCGAAGACATTTCCAGAGACTAATGACAAGCTTTTTCCA
CTGAGAAAATAAGAGGCAGGACGTTACCCAATGTCTGCATTGCTCAGTAGGCAGTTTTACACACAAAT
GCAGCACTGGAAGAGGTGCCAGATGCGTGAATGATATGCATAACTTTACTAACGAAATACTGACTTGG
GATCTGGTCCCCAAGTAAAACAACAGACCGATAAATCTTCTGACTGCATCACAAAACTACAGGACTT
GACATGAGTGAATATAATTCAGAAATCCAGTTGTAATCTTAAAAGTACTCATCAGAAAATCCT
GTATGTTCTATTGATGGGAGCACTCCCATCACTAAATCAACAGGTGATTGGGCAGAAAGCATCTCCAG
CAAGAAAATGCTATTACTGGAAAACCTGTACCTGACTCTTTGAATTCCACAAAAGAAATTCAGTATCAAA
GGTGTGCAAGGCAATATGCAGAAAACAGGATGACACACTTGCAAAATACCTGGATCTCCACCTGAGAAA
TGTGACTCTTTGGGTTCTGGAGTGGCCACAGTAAAAGTGGTTTTACCTGATGACCACCTGAAAGATGAA
ATGGACTGGCAGAGCTCTGCATTGGGAGAAAATCACAGAAGCTGATAGTTCTGGTGAAGTCTGATGACACA
GTAATAGAGGACATCACAGCAGATACATCATTTGAAAATAACAAAATTCAGGCTGAAAACCTGTTTCC
ATTCGAAGTCTGTTGTAAAACAGGTGAAAGAGAAAATCAAAGAGATCCAGTTGTGAGAGAGAAGAA
AAAACATCAAAAATTTGAAGAATTTGGTCAAGTACTGAGTGCATCAAGATCAGCTGATTTCTT
GGAAGGATCCAGCTAGTGAGGCAGCATGTTCAAAGTACCCGATACGAATGTCTCCTTAGAAGATGTG
AGTGAAGTTGCTCCTGAAAAGCCTATTACTACTGAGAACCCAAAATTCCTTCAACAGTGTCTCCAAT
GTTTTTAATGAGACAGAATTTCTCATTAAATGTGACAACATCTGCCTATTTGGAGTCATTACATGGGAAA
AATGTTAAACATATAGATGATTCCTCCCGAGGACCTGATAGCAGCCTTACAGAAACCAGAGATAAA
GGAATAGTAGATAGTAAAAGAAATGCTTTTAAAGCAATATCAGAGAAGATGACAGACTTTAAAACAAT
CCTCCTGTAGAAGTCTTACATGAAAATGAGTCCGGTGGTCTGAAATTAAGACATTGGAAGCAAATAC
AGTGAACAAAGCAAAGAAAACAAATGGAAGTGAAGCTCTAGGTGTTTTCCCTACCCAAGTACTCCAGTA
GCATCTTTGACTTAGAACAAGAACAGCTCACAATTAAGGCTCTTAAAGAATTAGGTGAAAGACAGGTT
GAGAAGTCAACTTCTGCACAGCGTGACGAGAATTGCCTTCTGAAGAAGTACTGAAGCAAATTTACA
TTTGCTCCAGAATCTTGCCACAGAGATCATATGACATCCTAGAACGTAATGTCAAGAATGGATCTGAT
CTTGGGATTTCCAGAAAGCCATCACTATCAGAGAACTACTAGGGTAGATGCTGTTTCCAGCCTTAGC
AAGACTGAATTGGTAAAAAGCATGTCCTAGCAAGACTTCTGACAGACTTCTCAGTGCACGATCTGATT
TTCTGGAGAGATGTGAAGAAGACTGGGTTTGTCTTTGGCACCACGCTGATCATGTGCTTTCCCTGGCA
GCTTTCAGTGTATCAGTGTGGTTTCTTACCTCATCCTGGCTCTTCTCTGTGACCATCAGCTTCAGG
ATCTACAAGTCCGTATCCAAGCTGTACAGAAGTCAGAAGAAGGCCATCCATTCAAAGCCTACCTGGAC
GTAGACATTACTCTGCCTCAGAAGCTTTCCATAATTACATGAATGCTGCCATGGTGCACATCAACAGG
GCCCTGAAAATCATTATTGCTCTTTCTGGTAGAAGATCTGGTTGACTCCTTGAAGCTGGCTGTCTTC
ATGTGGCTGATGACCTATGTTGGTGTGTTTTAACGGAATCACCCCTTCTAATTCTGCTGAACTGCTC
ATTTTCAGTGTCCGATTGTCTATGAGAAGTACAAGACCCAGATTGATCACTATGTTGGCATCGCCCGA
GATCAGACCAAGTCAATTGTTGAAAAGATCCAAGCAAATCCCTGGAATCGCCAAAAAAGGCAGAA
TAA
  
```

Restriction Sites: SgfI-MluI
ACCN: NM_001265590
Insert Size: 2763 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001265590.1
RefSeq Size:	4640 bp
RefSeq ORF:	2763 bp
Locus ID:	10313
UniProt ID:	O95197
Cytogenetics:	11q13.1
Protein Families:	Transmembrane
MW:	100.8 kDa
Gene Summary:	<p>This gene belongs to the reticulon family of highly conserved genes that are preferentially expressed in neuroendocrine tissues. This family of proteins interact with, and modulate the activity of beta-amyloid converting enzyme 1 (BACE1), and the production of amyloid-beta. An increase in the expression of any reticulon protein substantially reduces the production of amyloid-beta, suggesting that reticulon proteins are negative modulators of BACE1 in cells. Alternatively spliced transcript variants encoding different isoforms have been found for this gene, and pseudogenes of this gene are located on chromosomes 4 and 12. [provided by RefSeq, May 2012]</p> <p>Transcript Variant: This variant (6) includes an alternate exon in the coding region, but maintains the reading frame, compared to variant 1. The encoded isoform (f) is longer than isoform a.</p>