

Product datasheet for **SC332828**

FAM153B (NM_001265615) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: FAM153B (NM_001265615) Human Untagged Clone
Tag: Tag Free
Symbol: FAM153B
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC332828 representing NM_001265615.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGTGGACAAAGACACAGAGAGAGACATTGAGATGAAACGGCAACTACGGCGACTACGGGAGCTCCAC
CTATACAGCACATGGAAGAAGTACCAAGAGGCGATGAAGACATCCTTGGGAGTTCCACAATGTGAGCGT
GACGAAGGCTCCTTGGGCAAGCCATTGTGTCCACCCGAGATACTCTCGGAGACGTTGCCAGGCTCTGTG
AAGAAAAGGGTATGCTTTCATCAGAAGATCATCTAGAGGAGTTTATAGCAGAACATCTCCCTGAAGCA
TCCAATCAGAGTCTCCTCACTGTTGCCATGCAGACACAGGCATCCAAACCAACGGTGACCTGGAAGAC
CTGGAGGAGCATGGGCCAGGGCAGACAGTCTCTGAGGAAGCCACAGAAGTTCACATGATGGAGGGGGAC
CCAGACACACTGGCCGAACCTTCTGATCAGGGATGTACTTCAGGAGCTGTCCAGTTACAACGCGCAGGAG
GAGGACCCAGAGGAGGTGAAGACATCCTTGGGAGTTCACAACGTGGTGACCTGGAAGACCTGGAGGAG
CATGTGCCAGGGCAGACAGTCTCTGAGGAAGCCACAGGGGTTACATGATGCAGGTGGACCCAGCCACG
CCGGCAAAGAGTGACCTGGAAGACCTGGAGGAGCATGTGCCAGGGCAGACAGTCTCTGAGGAAGCCACA
GGGGTTTACATGATGCAGGTGGACCCAGCCACACTGGCAAAGCAATTGGAAGACTCCACCATTACAGGC
AGCCACCAGCAGATGTCAGCAAGTCCTTCTCTGCACCTGCAGAAGAAGCAACAGAAAAGACCAAAGTG
GAAGAGGAAGTAAAACAGAAAGCCCAAGAAGAAAACAGGAAGCCAGCAAGAAAAGCCGGTGAAT
GTCCTGAAATGTTGGGACATTTTAATATATT TAG
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001265615

Insert Size: 933 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).


[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001265615.1](#)

RefSeq Size: 1794 bp

RefSeq ORF: 933 bp

Locus ID: 202134

Cytogenetics: 5q35.2

MW: 34.7 kDa