

Product datasheet for **SC332846**

PLEKHB2 (NM_001267064) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: PLEKHB2 (NM_001267064) Human Untagged Clone
Tag: Tag Free
Symbol: PLEKHB2
Synonyms: EVT2
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC332846 representing NM_001267064.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
ATGGCGTTTGTGAAGAGTGGCTGGTTGCTGCGACAGAGTACTATTTTGAAGCGCTGGAAGAAGAACTGG
TTTGATCTGTGGTCGGATGGTCACCTGATCTATTATGATGACCAGACTCGGCAGAATATCGAGGATAAG
GTCCACATGCCAATGGACTGCATCAACATCCGCACGGGGCAGGAATGTCGGGATACTCAGCCCCGGAT
GGAAAGTCAAAGACTGCATGCTCCAGATTGTTTGTGCGAGATGGGAAAACAATTAGTCTTTGTGCAGAA
AGCACAGATGATTGCTTGGCCTGGAATTTACTCTCAAGATTCTAGGACAAACACAGCGTATGTGGGC
TCTGCAGTCATGACCGATGAGACATCCGTGGTTTCCTCACCTCCACCATACACGGCCTATGCTGCACCG
GCCCTGAGGTAGGGAGAACCCTGAGCCTCCAGGCTTATGGCTATGGGCCATACGGTGGTGCATACCCG
CCAGGAACTCAAGTTGTCTACGCTGCGAATGGGCAGGCGTATGCCGTGCCCTACCAGTACCCATATGCA
GGACTTTATGGACAGCAGCCTGTAACCAAGTCATCATTGCGAGAGCGCTATCGAGACAACGACAGCGAC
CTGGCACTGGGCATGCTGGCAGGAGCAGCCACGGGCATGGCCTTAGGGTCTCTATTTTGGGTCTTCTAG
```

Restriction Sites: Sgfl-Mlul
ACCN: NM_001267064
Insert Size: 690 bp
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001267064.1](#)

RefSeq Size: 4347 bp

RefSeq ORF: 690 bp

Locus ID: 55041

UniProt ID: [Q96CS7](#)

Cytogenetics: 2q21.1

MW: 25.5 kDa

Gene Summary: Involved in retrograde transport of recycling endosomes.[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (6) has multiple differences in the 3' coding region, compared to variant 4, which results in a protein (isoform 6) with a shorter C-terminus, compared to isoform 4.