

Product datasheet for **SC332865**

ZNF248 (NM_001267597) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: ZNF248 (NM_001267597) Human Untagged Clone
Tag: Tag Free
Symbol: ZNF248
Synonyms: bA162G10.3
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC332865 representing NM_001267597.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGAACAAATCCCAGGAACAAGTGTCAATCAAGGATGTATGTGTGGACTTCACTCAGGAAGAGTGGTAT
CTGCTGGACCCTGCTCAGAAGATTCTATACAGAGATGTGATCCTGGAAAATTATAGCAATCTTGTCTCA
GTAGGGTATTGCATTACTAAACCAGAAGTGATCTTTAAGATCGAGCAAGGAGAAGAGCCCTGGATATTA
GAAAAAGGATTCCAAGCCAGTGCCACCCAGAAAGGAAATGGAAAGTTGATGACGTGTTAGAGAGCAGC
CAGGAAAATGAAGATGACCATTTTTGGGAGCTTCTATTCCACAACAACAAAACAGTAAGTGTAGAAAAT
GGAGATAGAGGAAGCAAAACTTTCAATTTGGGCACAGACCCTGTTTCTTTAAGAAATTATCCCTATAAA
ATATGTGACTCATGTGAATGAATTTGAAAAATATTTCCGGCTTAATTATTAGTAAAAAGAACTGTTCC
AGAAAGAAGCCTGATGAGTTAATGTATGTGAGAAATTGCTCCTTGATATTAGGCATGAGAAAATCCCT
ATTGGAGAGAAGTCTTATAAATATGATCAAAAAAGGAATGCCATTAATTATCACCAGGATCTCAGTCAG
CCAAGTTTTGGCCAATCTTTTGGTATAGTAAAAATGGACAAGGCTTCCATGATGAGGCAGCATTTTTT
ACAAATAAGAGATCTCAGATAGGAGAGACAGTCTGTAATATAACGAATGTGGAAGAACCCTTATTGAA
AGTTTAAAGCTGAATATATCTCAAAGACCTCATTTGAAATGGAGCCGTATGGATGCAGTATTTGCGGG
AAGTCTTCTGCATGAATTTAAGGTTTGGACATCAGAGAGCTCTTACAAAGGACAATCCTTATGAATAT
AATGAATATGGGAAAATCTTCTGTGACAATTCAGCTTTTATTATCCATCAGGGAGCTTACACAAGAAAG
ATTCTCCGTGAATATAAAGTGAGTGACAAAACCTGGGAAAAGTCAGCTCTTTAAAACATCAAATAGTA
CACATGGGGGAAAAGTCTTATGATTACAATGAAAATGGGAGTAATTTAGCAAGAAGTCACATCTTACC
CAGCTTCGGAGAGCTCACACAGGAGAAAAAACCTTTGAATGTGGTGAATGTGGGAAAACCTTCTGGGAG
AAGTCAAACCTCACTCAACATCAGAGAACACACACAGGAGAGAAGCCCTATGAATGTACTGAATGTGGG
AAAGCCTTTTCCAGAAACCACCTGACCAACCATCAGCGAACACATACAGGAGAAAAACCTATGAA
TGTAAGCAATGTGAAAAACATTCTGTGTGAAGTCAAACCTCACTGAACATCAGAGAACACACACAGGG
GAGAAGCCCTATGAATGTAATGCATGTGGGAAATCCTTCTGCCACAGATCAGCCCTCACTGTGCATCAG
AGAACACACACAGGGGAGAAACCGTTTATATGTAATGAATGTGAAAAATCCTTCTGTGTGAAGTCAAAC
CTATTGTACATCAAAGAACTCACACTGGGAGAAACCATATAAGTGAATGAATGTGGGAAAACCTTC
TGTGAAAAATCAGCTCTCACTAAACATCAGAGGACTCACACAGGGGAGAAAGCCGTATGAGTGAATGCA
TGTGGGAAAGACCTTTAGTCAGAGGTCAGTGCTACCAAACATCAGAGAATTACACAAGGGTGAAAGCT
CTTTCAACATCCTGA
  
```

Restriction Sites: Sgfl-Mlul



[View online »](#)

ACCN:	NM_001267597
Insert Size:	1740 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001267597.1</u>
RefSeq Size:	4819 bp
RefSeq ORF:	1740 bp
Locus ID:	57209
UniProt ID:	<u>Q8NDW4</u>
Cytogenetics:	10p11.21
Protein Families:	Transcription Factors
MW:	67.1 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (2) differs in the 5' UTR, compared to variant 1. Variants 1, 2, and 6-11 all encode the same isoform (1).