

Product datasheet for **SC332869**

ARPP21 (NM_001267617) Human Untagged Clone

Product data:

| | |
|---------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | ARPP21 (NM_001267617) Human Untagged Clone |
| Tag: | Tag Free |
| Symbol: | ARPP21 |
| Synonyms: | ARPP-21; R3HDM3; RCS; TARPP |
| Vector: | pCMV6-Entry (PS100001) |



[View online »](#)

Fully Sequenced ORF: >SC332869 representing NM_001267617.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGTCTGAGCAAGGAGACCTGAATCAGGCAATAGCAGAGGAAGGAGGGACTGAGCAGGAGACGGCCACT
CCAGAGAACGGCATTGTTAAATCAGAAAGTCTGGATGAAGAGGAGAACTGGAAGTGCAGAGGCGGCTG
GAGGCTCAGAATCAAGAAAGAAGAAAATCCAAGTCAGGAGCAGGAAAAGTAACTGACTCGCAGCCTT
GCTGTCTGTGAGGAATCTTCTGCCAGACCAGGAGGTGAAAGTCTTCAGGATCAGGAATCAATTCATTTA
CAGCTTTCCAGTTTTTCCAGCCTGCAAGAGGAGGATAAATCTAGGAAAGATGACTCTGAAAGAGAAAAA
GAAAAGGATAAAAAACAAAGATAAAAACCTCTGAAAAACCAAGATCAGAATGTTATCAAAAAGATTGCAGC
CAAGAATACACGGATTCTACAGGCATAGACTTACACGAGTTTCTGATTAACACATTAAGAATAATTCC
AGGGACAGGATGATACTTTTAAAATGGAGCAGGAAATTATTGATTTTATTGCTGACAACAATAATCAT
TATAAAAAGTTCCCTCAGATGTCATCGTATCAGAGGATGCTTGTCCATCGAGTGGCAGCTTATTTTGA
TTGGATCACAATGTGGATCAAACAGGAAAATCTGTTATCATCAACAAGACCAGCAGCACCAGAATACCA
GAGCAAAGGTTTTGTGAACATTTAAAAGATGAAAAGGTGAAGAATCCCAGAAGCGGTTTATCTTGAAG
CGAGATAACTCTAGTATTGATAAAGAAGACAATCAGTCAGTTTGCTCCCAGGAAAGCCTTTTTGTGGAA
AACAGGGGCAACAGAGATGGCTCAGGGAGAACATCTGGGAGTCGACAGAGCAGCTCAGAAAATGAACTC
AAGTGGTCTGACCACAAAGGGCCTGGAGCAGCACAGACTCCGACAGTCCAACCGCAATCTAAAGCCC
GCCATGACCAAGACGGCAGTTTTGGGGGCATCACGGTGTGACCAGGGGTGACAGCACTTCCAGTACT
AGGAGTACCGGAAGCTGTCCAAGCAGGTTCCGAGTCTTCCAGCAGTGCAGGCTCCTCAGGATCGCTG
TCCCGCACCCATCCACCTCTCCAGAGCACACCCCTAGTCTCAGGTGTGGCAGCTGGCTCTCCAGGCTGT
GTGCCTTATCCAGAGAAATGGAATAGGGGGCCAGGTTGCTCCCAGCAGCACCAGCTACATCCTCCTCCA
ATTCCCAGTGAACAGGCATCCCGCCTGGAAGCATCCTTCTTAATCCACACACAGGCCAGCCCTTTGTG
GGGCAGTCCCAACAGCAGCCACCACAGCAGCAGCCCTCCCGCAGCCCAACAGCAGGTCCAGCCACCG
CAGCCACAGATGGCAGGCCCTCTGGTCACTCAGTCTGTCCAGGGGCTGCAGGCTTCTCCAGTCACTG
CAATATCCAGCAGTCTCTTTTCTCCCGCAGCCTCCTACCTGTGTCTCCAACGCAGCACTTTCCCATG
AGAGATGATGTGGCAACACAGTTTGGCCAGATGACCTGAGCCGGCAGTCTCGGGGGAGACTCCTGAA
CCCCCATCAGGTCTGTCTACCCATCCTCCCTTATGCCACAGCCGGCCAGCAGCCAGCTATGTAATC
GCCTCTACAGGCCAGCAGCTTCTACAGGAGGATTCTCAGGCTCTGGCCCTCCCCTCCTCCAGCAGGTC
CTCCAGCCCCCTCCCTACCACAGGGATTTGTGCAACAGCCTCCGCTGCACAGATGCCTGTATATTAT
TACCCATCTGGTCACTACCTCAACCACGCAACAGTACCGCCCATGGCCCGGTTAGTACAAC
GCTCAGAGGAGTCAACAGATGCCACAGGAGCAGCAAGCAAGCAGGTTACCAGCCAGTCTTGTCTGGTCAA
CAGGGATTCCAAGGCCTAATAGGAGTGCAGCAGCCACCTCAGAGTCAGAACGTGATAAATAACCAACAA
GGAATCCGGTGCAAAGCGTGATGGTTTTCTACCCAACAATGTCTTCTTATCAGGTGCCAATGACCCAG
GGTTCTCAAGGACTGCCCCAGCAGTCAACCAACAGCCAATCATGCTACCTAACCAGGCAGGTCAAGGG
TCACTCCCAGCCACTGGAATGCCTGTTACTGTAATGTACACCGCCACCCCTCAGAACAACCTTAGG
CTGATTGGCCCACTGCCCTCCAGCACTGTCCAGTGTGACAGTGTGACAGTGTGACAGTGTGACAGTGTG
AGTATGAGCAATGCTGGTTGGCAGGTCAAATTC
```

Restriction Sites: SgfI-MluI

ACCN: NM_001267617

Insert Size: 2382 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001267617.1](#)

RefSeq Size: 3184 bp

RefSeq ORF: 2382 bp

Locus ID: 10777

UniProt ID: [Q9UBL0](#)

Cytogenetics: 3p22.3

MW: 86.1 kDa

Gene Summary: This gene encodes a cAMP-regulated phosphoprotein. The encoded protein is enriched in the caudate nucleus and cerebellar cortex. A similar protein in mouse may be involved in regulating the effects of dopamine in the basal ganglia. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Jun 2012]
Transcript Variant: This variant (6) contains a distinct 5' UTR and lacks an in-frame exon in the central coding region, compared to variant 1. The resulting isoform (3) is shorter, compared to isoform 1.