

Product datasheet for **SC332887**

CLOCK (NM_001267843) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CLOCK (NM_001267843) Human Untagged Clone
Tag:	Tag Free
Symbol:	CLOCK
Synonyms:	bHLHe8; KAT13D
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC332887 representing NM_001267843.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGTTGTTTACCCTAAGCTGTAGTAAAAAGAGCTCGATTGTTGACAGAGATGACAGTAGTATTTTGTAT
GGTTGGTGGAAAGATGACAAGGACAAAGCGAAAAGAGTATCTAGAAACAACTGAAAAAGAACGT
AGAGATCAATTTAATGTTCTCATTAAAGAACTGGGATCCATGCTTCTGTGTAATGCTAGAAAGATGGAC
AAATCTACTGTTCTGCAGAAAAGCATTGATTTTTTACGAAAACATAAAGAAATCACTGCACAGTCAGAT
GCTAGTGAATTCGACAGGACTGAAAACCTACATTCCTTAGTAATGAAGAGTTTACACAATTAATGTTA
GAGGCTCTTGATGGTTTTTTTTAGCAATCATGACAGATGGAAGCATAATATATGTGTCTGAGAGTGTA
ACTTCATTACTTGAACATTTACCATCTGATCTTGTGGATCAAAGTATATTTAATTTTATCCAGAAGGG
GAACATTCAGAGGTTTAAAAACTCTCTACTCATCTGCTGGAAAGTGATTCATTAACCCAGAATAT
TAAAAATCAAAAATCAGTTAGAATTCTGTTGTCACATGCTGCGAGGAACAATAGACCCAAAGGAGCCA
TCTACCTATGAATATGTAATTTATAGGAAATTTCAAATCTTTAACAGTGTATCCTCTTCAGCACAC
AATGGTTTTGAAGGAACTATACAACGCACACATAGGCCATCTTATGAAGATAGAGTTTGTGTTGAGCT
ACTGTCAGGTTAGCTACACCTCAGTTCATCAAGGAAATGTGCACTGTTGAAGAACCAATGAAGAGTTT
ACATCTAGACATAGTTTGAATGGAAGTTTCTGTTTCTAGATCACAGGGCACCACCATAATAGGGTAT
TTGCCATTTGAAGTTCTGGGAACATCAGGCTATGATTACTATCATGTGGATGACCTAGAAAAATTTGGCA
AAATGTCATGAGCACTTAATGCAATATGGGAAAGGCAAATCATGTTATTATAGGTTCTGACTAAGGGG
CAACAGTGGATTTGGCTCAGACTCATTATTATCACTTACCATCAGTGAATTAAGGCCAGAGTTT
ATTGTTGTACTCACACTGTAGTAAGTTATGCAGAAGTTAGGGCTGAAAGACGACGAGAAGTTGGCATT
GAAGAGCTCTTCTGAGACAGTGTGACAAAAGCCAAGATTCTGGGTCAGATAATCGATAAACACA
GTCAGTCTCAAGGAAGCATTGGAAGGTTTGATCACAGCCCAACCCCTTCTGCCTTCTCGGAGTTCA
AGAAAATCATCTCACACGGCCGCTCAGACCCTTCTCAACACCAACCAAGATCCCGACGGATACGAGC
ACTCCACCAGGCAGCATTACCAGCTCATGAGAAGATGGTGCAAGAAGGTCATCATTTAGTAGTCAG
TCCATAAATTCAGTCTGTTGGTTTCATTAACACAGCCAGTGTGCTCAAGCTACAAATTTACCA
ATTCCACAAGGCATGTCCAGTTTCAGTTTTTCAGCTCAATTAGGAGCCATGCAACATCTGAAAGACCAA
TTGGAACAACGGACACGCATGATAGAAGCAAATTTTCATCGGCAACAAGAAGAACTAAGAAAAATCAA
GAACAACCTCAGATGGTCCATGGTCAGGGGCTGCAGATGTTTTTGCAACAATCAAATCCTGGGTTGAAT
TTTGGTTCGGTTCAACTTCTTCTGAAAATTCATCTAATATCCAGCAACTGCACCTATAAATATGCAA
GGCCAAGTTGTTCTACTAACCAGATTCAAAGTGGAAATGAATACTGGACACATTGGCACAACCTCAGCAC
ATGATACAACAACAGACTTTACAGAGTACATCAACTCAGAGTCAACAAAATGACTGAGTGGGCACAGT
CAGCAAACATCTCTACCCAGTCAGACACAGAGCACTTTACAGCCCACTGTATAACACTATGGTGATT
TCTCAGCTGCAGCCGGAAGCATGGTCCAGATTCATCTAGTAGCCACAAAACAGCACCCAGAGTGCT
GCAGTAACTACATTCACTCAGGACAGGCAGATAAGATTTTCTCAAGGTCAACAACCTGTGACCAAATTA
GTGACTGCTCCTGTAGCTTGTGGGCAGTCATGGTACCTAGTACTATGCTTATGGCCAGGTGGTGACT
GCATATCTACTTTTGTACACAACAGCAACAGTCACAGACATTGTCAGTAACGCAGCAGCAGCAGCAG
CAGAGTCCCAGGAGCAGCAGCTCACTTCAGTTCAGCAACCATCTCAGGCTCAGCTGACCCAGCCACCG
CAACAATTTTACAGACTTCTAGGTTGCTCCATGGGAATCCCTCAACTCAACTCATTCTCTCTGTGCA
TTTCTCTACAACAGAGCACCTTCCCTCAGTCACATCACAGCAACATCAGTCTCAGCAACAGCAGCAA
CTCAGCCGGCACAGGACTGACAGCTTGCCCGACCCCTTCAAGGTTCAACCCAGTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-Mlul
ACCN: NM_001267843
Insert Size: 2541 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001267843.1](#)

RefSeq Size: 10981 bp

RefSeq ORF: 2541 bp

Locus ID: 9575

UniProt ID: [O15516](#)

Cytogenetics: 4q12

Protein Families: Druggable Genome, Transcription Factors

Protein Pathways: Circadian rhythm - mammal

MW: 95.3 kDa

Gene Summary:

The protein encoded by this gene plays a central role in the regulation of circadian rhythms. The protein encodes a transcription factor of the basic helix-loop-helix (bHLH) family and contains DNA binding histone acetyltransferase activity. The encoded protein forms a heterodimer with ARNTL (BMAL1) that binds E-box enhancer elements upstream of Period (PER1, PER2, PER3) and Cryptochrome (CRY1, CRY2) genes and activates transcription of these genes. PER and CRY proteins heterodimerize and repress their own transcription by interacting in a feedback loop with CLOCK/ARNTL complexes. Polymorphisms in this gene may be associated with behavioral changes in certain populations and with obesity and metabolic syndrome. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2014]

Transcript Variant: This variant (1) represents the longer transcript. Both variants 1 and 2 encode the same protein. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.