

Product datasheet for **SC332935**

RTBDN (NM_001270443) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: RTBDN (NM_001270443) Human Untagged Clone
Tag: Tag Free
Symbol: RTBDN
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC332935 representing NM_001270443.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
ATGGACTGCAGGGTCCACATGCGACCCATCGGCCTGACGTGGGTGCTGCAACTGACCTTGGCATGGATC
CTGCTAGAAGCCTGTGGAGGGAGCCGCCACTCCAAGCCAGGTCCCAGCAACACCATGGGCTGGCAGCT
GATCTGGGCAAAGGCAAGCTGCACCTGGCAGGACCTTGTGTCCCTCAGAGATGGACACAACAGAGACA
TCGGGCCCTGGAAACCATCCAGAACGCTGTGGAGTGCCGAGCCCTGAATGCGAATCCTTCCTGGAACAC
CTCCAACGTGCCCTTCGCAGTCGCTTCCGCCTGCGGCTATTGGGGTACGCCAGGCACAGCCGCTCTGC
GAGGAGCTCTGCCAGGCCTGGTTCGCCAACTGCGAAGATGATATCACCTGCGGCCGACTTGGCTCCCA
CTCTCAGAAAAAAGGGGCTGTGAGCCCAGCTGCCTTACCTATGGACAGACCTTCGCAGACGGGACGGAC
CTTTGTGCTCGGCTCTGGGCCACGCCCTACCGGTGGCTGCTCCTGGAGCCCGTCACTGCTTCAACATC
TCCATCTCCGCGGTACCTCGTCCAGACCAGGACGACGGGGCCGGGAAGCTCCCTCCCGCGCTTCCCGC
AGCCCTCGCACCTCCATCCTGGACGCTGCGGGCAGCGGGAGTGGCAGTGAAGCGGCAGCGGCCCTAG
```

Restriction Sites: SgfI-MluI
ACCN: NM_001270443
Insert Size: 690 bp
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001270443.1](#)

RefSeq Size: 1018 bp

RefSeq ORF: 690 bp

Locus ID: 83546

UniProt ID: [Q9BSG5](#)

Cytogenetics: 19p13.13

Protein Families: Druggable Genome, Secreted Protein

MW: 24.6 kDa

Gene Summary: This gene was first identified in a study of human eye tissues. The protein encoded by this gene is preferentially expressed in the retina and may play a role in binding retinoids and other carotenoids as it shares homology with riboflavin binding proteins. Alternative splicing results in multiple transcript variants and protein isoforms. [provided by RefSeq, Jul 2012]
Transcript Variant: This variant (6) differs in the 5' UTR, compared to variant 1. Variants 1, 6 and 7 encode the same isoform (1).