

Product datasheet for **SC332941**

UBE3B (NM_001270450) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: UBE3B (NM_001270450) Human Untagged Clone
Tag: Tag Free
Symbol: UBE3B
Synonyms: BPIDS; KOS
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC332941 representing NM_001270450.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
ATGTTTACCCTGTCTCAGACCTCGAGAGCATGGTTCATCGATAGAGCCCGTCAGGCACGAGAAGAAAGG
CTTGTGTCAGAAGGAACGGGAGCGGGCAGCTGTTGTGATCCAGGCCATGTCCGGAGTTTTCTGTCCGG
AGTCGACTGCAGAGAGATATCAGGAGAGAGATTGATGACTTTTTAAAGCAGATGACCCTGAGTCCACT
AAAAGAAGTGCACCTTTGATTTTCAAGATTGCCAGGAACTGCTGTTCTATTCAAGATCAAGAGGAT
AATGAGAGATTTGAGAAGTTGTGTCGCAGCATCCTGAGCAGCATGGATGCTGAGAATGAGCCTAAGGTG
TGGTATGTGTCCTGGCTTGTCTAAGGACCTCACCTCCTTTGGATTCAACAGATCAAGAACATTTTG
TGGTACTGCTGTGATTTTCTCAAGCAGCTCAAGCCTGAAATCCTGCAGGACTCCCGACTCATCACCTG
TACCTCACGATGCTTGTACCTTACAGACACTTCAACGTGAAAATTCTTCGGGAAAAGGTGAAAGT
CTTCGACCAGCGATGAACCACATTTGTGCAAAATAATGGGACATCTCAACCAGCATGGATTTATTCT
GTGCTGCAGTGTGATGGGCTGTTTCTGATTTGGTTTCATATGCTCCTCACAAACCTGTGAGG
TGGTCCGTTGGCAGAAGCTGGTATGACTGGCAGTTGTCTCGCTAG
```

Restriction Sites: Sgfl-MluI
ACCN: NM_001270450
Insert Size: 735 bp
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001270450.1](#)

RefSeq Size: 1146 bp

RefSeq ORF: 735 bp

Locus ID: 89910

UniProt ID: [Q7Z3V4](#)

Cytogenetics: 12q24.11

Protein Families: Druggable Genome

Protein Pathways: Ubiquitin mediated proteolysis

MW: 28.9 kDa

Gene Summary: The modification of proteins with ubiquitin is an important cellular mechanism for targeting abnormal or short-lived proteins for degradation. Ubiquitination involves at least three classes of enzymes: E1 ubiquitin-activating enzymes, E2 ubiquitin-conjugating enzymes, and E3 ubiquitin-protein ligases. This gene encodes a member of the E3 ubiquitin-conjugating enzyme family which accepts ubiquitin from an E2 ubiquitin-conjugating enzyme and transfers the ubiquitin to the targeted substrates. A HECT (homology to E6-AP C-terminus) domain in the C-terminus of the longer isoform of this protein is the catalytic site of ubiquitin transfer and forms a complex with E2 conjugases. Shorter isoforms of this protein which lack the C-terminal HECT domain are therefore unlikely to bind E2 enzymes. Alternatively spliced transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, Jul 2012]

Transcript Variant: This variant (5) uses an alternate splice site in the 5' UTR and contains an alternate 3' exon, compared to variant 1. These differences result in a distinct 3' UTR and 3' coding region and a protein (isoform 3) with a shorter and distinct C-terminus, compared to isoform 1. Variants 4-6 encode the same protein (isoform 3).