

Product datasheet for **SC332945**

WWP2 (NM_001270454) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	WWP2 (NM_001270454) Human Untagged Clone
Tag:	Tag Free
Symbol:	WWP2
Synonyms:	AIP2; WWp2-like
Vector:	pCMV6-Entry (PS100001)



[View online »](#)

Fully Sequenced ORF: >SC332945 representing NM_001270454.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGCATCTGCCAGCTCTAGCCGGGCAGGAGTGGCCCTGCCTTTTGAGAAGTCTCAGCTCACTTTGAAA
GTGGTGTCCGCAAAGCCCAAGGTGCATAATCGTCAACCTCGAATTAACCTCCTACGTGGAGGTGGCGGTG
GATGGACTCCCCAGTGAGACCAAGAAGACTGGGAAGCGCATTGGGAGCTCTGAGCTTCTCTGGAATGAG
ATCATCATTTTTGAATGTCACGGCACAGAGTCATTTAGATTTAAAGTCTGGAGCTGCCATACCTTGAGA
AATGAACCTGCTAGGCACCGCATCTGTCAACCTCTCCAACGTCTTGAAGAACAATGGGGGCAAAATGGAG
AACATGCAGCTGACCCTGAACCTGCAGACGGAGAACAAGGCAGCGTTGTCTCAGCGGAGAGCTGACA
ATTTTCTGGACGGGCCAACTGTTGATCTGGGAAATGTGCCTAATGGCAGTGCCTGACAGATGGATCA
CAGCTGCCTTCGAGAGACTCCAGTGAACAGCAGTAGCTCCAGAGAACCAGCAGCCCCCAGCACA
AACTGCTTTGGTGAAGATCCCGGACGCACAGACATTCGGGTGCTTCAGCCAGAACAACCCAGCAACC
GGCGAGCAAAGCCCCGGTCTCGGAGCCGGCACCAGCCGTCAGAAGTCAAGCCACAGTGGCTTG
GCCAATGGCACAGTGAATGATGAACCCACAACAGCCACTGATCCCAGAAGAACCTCCGTTGTTGGTGTG
ACGTCCCACCTGCTGCACCCTTGAGTGTGACCCCGAATCCCAACACGACTTCTCTCCCTGCCCCAGCC
ACACCCGGCTGAAGGAGAGGAACCCAGCACTTCGGGTACACAGCAGCTCCCAGCGGCTGCCAGGCCCCC
GACGCTCTGCCTGCTGGATGGGAACAGCGAGAGCTGCCAACGGACGTGTCTATTATGTTGACCACAAT
ACCAAGACCACCCTGGGAGCGGCCCTTCTCCAGGCTGGGAAAAACGCACAGATCCCCGAGGCGAGG
TTTTACTATGTGGATCACAATACTCGGACCACCACCTGGCAGCGTCCGACCGCGGAGTACGTGCGCAAC
TATGAGCAGTGGCAGTCGACGGGAATCAGCTCCAGGGGGCCATGCAGCACTTCAGCCAAAGATTCCTC
TACCAGTCTTCGAGTCTTCGACTGACCATGATCCCCTGGGCCCTCCCTCCTGGCTGGGAGAAGAGA
CAGGACAATGGACGGGTGATTACGTGAACCATAACACTCGCACGACCCAGTGGGAGGATCCCCGGACC
CAGGGGATGATCCAGGAACCGCTCTGCCCCAGGATGGGAGATGAAATACACCAGCGAGGGGTGCGA
TACTTTGTGGACCACAATACCCGACCCACCTTTAAGGATCCTCGCCCGGGGTTTGAGTCGGGGACG
AAGCAAGTTCCCTGCTTATGACCGCAGTTTTCGGTGGAAGTATCACCAAGTCCGTTTCTCTGC
CATTCAAATGCCCTACCTAGCCACGTGAAGATCAGCGTTTCCAGGCAGACGCTTTTTCGAAGATTCCTTC
CAACAGATCATGAACATGAAACCCTATGACCTGCGCCCGCGGCTCTACATCATATGCGTGGCGAGGAG
GGCCTGGACTATGGGGGCATCGCCAGAGAGTGGTTTTTCTCCTGTCTCATGAGGTGCTCAACCCTATG
TATTGTTTATTTGAATATGCCGGAAGAACAATTACTGCCTGCAGATCAACCCCGCTCCTCCATCAAC
CCGGACCACCTCACCTACTTTTCGCTTTATAGGCAGATTCATCGCCATGGCGTGTACCATGGAAAGTTC
ATCGACACGGGCTCACCTCCCTTTTACAAGCGGATGCTCAATAAGAGACCAACCCTGAAAGACCTG
GAGTCCATTGACCTGAGTCTACAACCTCATTGTCTGGATCAAAGAGAACAACCTGGAAGAATGTGGC
CTGGAGCTGTACTTCATCCAGGACATGGAGATACTGGGCAAGGTGACGACCCACGAGCTGAAGGAGGGC
GGCGAGAGCATCCGGGTACAGAGGAGAACAAGGAAGGTACATCATGCTGCTGACTGACTGGCGTTTC
ACCCGAGGCGTGAAGAGCAGACCAAAGCTTCTGGAATGGCTTCAACGAGGTGGCCCCGCTGGAGTGG
CTGCGCTACTTTGACGAGAAAGAGCTGGAGCTGATGCTGTGCGGCATGCAGGAGATAGACATGAGCGAC
TGGCAGAAGAGCACCATCTACCGCACTACACCAAGAACAGCAAGCAGATCCAGTGGTCTGGCAGGTG
GTGAAGGAGATGGACAACGAGAAGAGGATCCGGCTGCTGCAGTTTGTACCGGTACCTGCCGCTGCC
GTCGGGGGATTTGCCGAACCTCATCGGTAGCAACGACACAGCAAGTTTTCATTGACAAAGTTGGCAAG
GAAACCTGGCTGCCAGAAGCCACACCTGCTTCAACCGTCTGGATCTTCCACCTACAAGAGCTACGAA
CAGCTGAGAGAGAAGCTGCTGTATGCCATTGAGGAGACCGAGGGCTTTGGACAGGAGTAA
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001270454

Insert Size: 2613 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001270454.1</u>
RefSeq Size:	4555 bp
RefSeq ORF:	2613 bp
Locus ID:	11060
UniProt ID:	<u>O00308</u>
Cytogenetics:	16q22.1
Protein Families:	Druggable Genome
Protein Pathways:	Ubiquitin mediated proteolysis
MW:	98.9 kDa
Gene Summary:	<p>This gene encodes a member of the Nedd4 family of E3 ligases, which play an important role in protein ubiquitination. The encoded protein contains four WW domains and may play a role in multiple processes including chondrogenesis and the regulation of oncogenic signaling pathways via interactions with Smad proteins and the tumor suppressor PTEN. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene, and a pseudogene of this gene is located on the long arm of chromosome 10. [provided by RefSeq, Jul 2012]</p> <p>Transcript Variant: This variant (5) differs in the 5' UTR, compared to variant 1. Variants 1 and 5 encode the same isoform (WWP2-FL).</p>