

Product datasheet for **SC333011**

JAKMIP2 (NM_001270941) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	JAKMIP2 (NM_001270941) Human Untagged Clone
Tag:	Tag Free
Symbol:	JAKMIP2
Synonyms:	JAMIP2; NECC1
Vector:	pCMV6-Entry (PS100001)



[View online »](#)

Fully Sequenced ORF: >SC333011 representing NM_001270941.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGTCCAAGAAAGGGCGAAATAAGGGCGAGAAGCCGAGGCACTCATTGTTGCCCTTCAAGCTGCCAAT
GAAGACCTCAGGACCAAGCTCACAGACATTCAGATAGAGCTGCATCAAGAGAAGTCCAAGGTATCAAAG
CTTGAAAGAGAGAAGACTCAAGAAGCGAAGAGGATTCGTGAGCTGGAGCAGCGCAAGCACACGGTGCTG
GTGACAGAACTCAAAGCCAAGCTCCATGAGGAGAAGATGAAGGAGCTGCAGGCTGTGAGGGAGAACCCTT
ATCAAGCAGCACGAGCAGGAAATGTCAAGGACGGTGAAGGTACGTGATGGAGAGATCCAGAGGCTCAAG
TCTGCTCTCTGTGCTCTCCGCGACGCGCAGCAGTGACAAAGTAAGGACAGCGCTCACCAATTGAGGCCCGG
GAGGAGGCCCGAAACTGTTTGACACAGAGCGCCTTAAGCTCTTACAGGAAATTGCGGACCTGAAAACG
GCCAAGAAGCAGGTGGACGAGGCTCTGAGCAATATGATCCAAGCAGATAAAATCAAGGCTGGGGACCTT
CGGAGTGAGCATCAGTCCCACCAAGAAGCCATCTCGAAGATCAAGTGGGAGTCGGAGCGGGATATTCGG
AGGCTGATGGATGAAATCAAAGCCAAGGACAGGATCATCTTTCCCTGGAAAAGGAACTGGAGACCCAG
ACAGGCTATGTACAGAACTCCAACCTCAGAAGGAGGCTTTGGACGAACAACTCTTCTGGTCAAGGAG
GCTGAGTGCAACATGAGCAGCCAAAACGAGAAATCCAGGAAGGCGAGGTGATGTTCCGAACACTGC
AGCAGTCTGATTTGCGAAGAAATCAAAAGAGAATAGCTGAATTGAATGCCACTATAAGAAAATTAGAA
GACAGGAATACCTTGCTTGAGATGAACGAAATGAACTGTTAAAACGTGTGCGGGAAACCAGAAAAGCAA
TGTAACCTCTCCTGGAAAGGAACAAGTGCCTCGCCAAGAGAAACGATGAACTGATGGTGTCTTGCAG
CGCATGGAAGAAAAACTAAAAGCCGTTACCAAGGAAAAATCAGAAATGAGAGAAAAATAACATCCCAT
CCACCCCTGAAGAAATTAATCTCTGAATGACCTCGACCAAGCTAATGAAGAACAAGAAACAGAGTTT
CTAAAACCTCAGGTCATTGAGCAACAGAACATTATTGATGAGCTCACAAGGGACCGAGAAAAGCTCATC
CGTAGAAGAAAGCATAGAAGAAGTTCCAAGCCAATTAAGAGGCCGTGTTTGGACCCGTTTATTGGCTAT
GATGAGGACTCTATGGATTCAGAGACATCATCCATGGCCTCATTTAGAACAGACAGAACCAGCTACT
CCTGATGATGACTTGGATGAAAGTTTAGCAGCTGAAGAATCTGAACTAAGATTTGACAATTAACAAAA
GAATATCAGGCCCTCAAAGAGCATATGCCCTCCTACAGGAGCAGACGGGAGGCATCATCGACGCTGAA
CGAGAAGCCAAGGCTCAAGAACAGCTCCAAGCAGAGGTGTAAGGTATAAAGCCAAAATTGAAGACCTG
GAAGCGACTCTGGCTCAGAAAGGCGAGGATTCACACTGGGTAGAAGATAAACAACTTTTTCATTAAGAGA
AACCAGGAGCTTTTAGAAAAGATAGAAAAACAGGAGGCAGAAAATCACCGGTTACAACAAGAACTACAG
GACGCCAGAGACCAGAATGAGCTGCTGGAGTTTCGAAACCTAGAGCTAGAAGAGAGAGAGACGATCC
CCTCCATTTAATCTCAAATTCACCCATCTCAGATGGTGTGAGTGTCTACAGATCTACTGTATGAAA
GAAGGTGTTAAGGATGTGAACATCCCTGATCTCATAAAGCAGCTTGATATCTTGGGTGATAATGGGAT
TTAAGAAATGAAGAACAAGTGGCCATAATTCAGGCCAGCACTGTGCTGTCCCTGGCAGAGAAGTGGATC
CAGCAGATTGAAGGAGCTGAGGCTGCCCTACACCAGAAAATGATGGAATTGGAAGTGACATGGAACAG
TTCTGCAAAAATAAAGGCTATCTGGAGGAAGAACTAGACTACAGAAAACAAGCTCTTGACCAAGCATAT
ATGAGAAATCCAGGAACTAGAAGCTACTTTGTACAATGCTCTACAGCAAGAAACTGTTATCAAGTTTGGT
GAATTATTAAGTGAAAAACAGCAAGAGGAGCTGAGGACGGCAGTAGAAAAGTTACGGCGGCAATGCTG
AGGAAGAGCAGAGAGTATGACTGTCAGATTCTTCAGGAGAGAATGGAGCTTTACAGCAAGCCCATCAG
AGAATTCGTGACTTAGAAGATAAAAACAGACATCCAGAAAAGACAAATAAAAGACTTAGAAGAAAAGTTT
CTGTTTCTATTCTGTTCTTCTCTTTGCCCTTATTCTATGGCCTTGA
  
```

Restriction Sites: Sgfl-Mlul

ACCN: NM_001270941

Insert Size: 2463 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001270941.1</u>
RefSeq Size:	9272 bp
RefSeq ORF:	2463 bp
Locus ID:	9832
UniProt ID:	<u>Q96AA8</u>
Cytogenetics:	5q32
MW:	96.3 kDa
Gene Summary:	<p>The protein encoded by this gene is reported to be a component of the Golgi matrix. It may act as a golgin protein by negatively regulating transit of secretory cargo and by acting as a structural scaffold of the Golgi. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2012]</p> <p>Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>