

Product datasheet for **SC333017**

IST1 (NM_001270975) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: IST1 (NM_001270975) Human Untagged Clone
Tag: Tag Free
Symbol: IST1
Synonyms: CHMP8; OLC1
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC333017 representing NM_001270975.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
ATGCTGGGCTCTGGATTTAAAGCTGAGCGCTTAAGAGTGAATTTGAGATTAGTCATAAATCGCCTTAAA
CTATTGGAGAAAAAGAAAACGGAACCTGGCCAGAAAGCAAGGAAGGAGATTGCTGACTATCTGGCTGCT
GGGAAAGATGAACGAGCTCGGATCCGTGTGGAGCACATTATCCGGGAAGACTACCTCGTGGAGGCCATG
GAGATCCTGGAGCTGACTGTGACCTGCTGCTGGCTCGGTTTGGCCTTATCCAGTCTATGAAGGAACCTA
GATTCTGGTCTGGCTGAATCTGTGTCTACATTGATCTGGGCTGCTCCTCGACTCCAGTCAGAAGTGGCT
GAGTTGAAAATAGTTGCTGATCAGCTCTGTGCCAAGTATAGCAAGGAATATGGCAAGCTATGTAGGACC
AACCAGATTGGAACGTGAATGACAGGCTAATGCACAAGCTGAGTGTGGAAGCCCAACCAAAATCCTG
GTGGAGAGATACCTGATTGAAATTGCAAAGAATTACAACGTACCCATGAACTGACTCTGTGGTCATG
GCAGAAGCTCCTCCTGGGGTAGAGACAGATCTTATTGATGTTGGATTACAGATGATGTGAAGAAAGGA
GGCCCTGGAAGAGGAGGGAGTGGTGGCTTACAGCACCAGTTGGTGGACCTGATGGAACGGTGCCAAATG
CCCATGCCATGCCATGCCTATGCCATCTGCAAATACGCCTTTCTCATATCCACTGCCAAAGGGACCA
TCAGATTTCAATGGACTGCCAATGGGGACTTATCAGGCCTTTCCAATATTCATCCACCTCAGATACCA
GCAACTCCCCATCGTATGAATCTGTAGATGACATTAATGCTGATAAGAATATCTCTTCTGCACAGATT
GTTGGTCTGGACCAAGCCAGAAGCCTCTGCAAAGCTTCCCTCCAGACCTGCAGATAACTATGACAAC
TTTGTCTACCAGAGTTGCCATCTGTGCCAGACACACTACCAACTGCATCTGCTGGTGCCAGCACCTCA
GCATCTGAAGACATTGACTTTGATGATCTTTCCGGAGGTTTGAAGAGCTGAAAAAGAAAACATAG
```

Restriction Sites: SgfI-MluI
ACCN: NM_001270975
Insert Size: 1101 bp
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001270975.1
RefSeq Size:	4039 bp
RefSeq ORF:	1101 bp
Locus ID:	9798
UniProt ID:	P53990
Cytogenetics:	16q22.2
MW:	40 kDa
Gene Summary:	<p>This gene encodes a protein with MIT-interacting motifs that interacts with components of endosomal sorting complexes required for transport (ESCRT). ESCRT functions in vesicle budding, such as that which occurs during membrane abscission in cytokinesis. There is a pseudogene for this gene on chromosome 19. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Aug 2012]</p> <p>Transcript Variant: This variant (1) encodes isoform a (also known as hIST1c). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>