

Product datasheet for **SC333057**

CCDC177 (NM_001271507) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CCDC177 (NM_001271507) Human Untagged Clone
Tag:	Tag Free
Symbol:	CCDC177
Synonyms:	C14orf162; PLPL
Vector:	pCMV6-Entry (PS100001)



[View online »](#)

Fully Sequenced ORF: >SC333057 representing NM_001271507.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGTGGACCCAGTGCCTGAAGAAGAGAAGGCAGGAGCGGAACCCGGCGACTCTGGAGGGGACGAGGCC
GTGGCGTCCGTGCCCCCTGATCCCAGGGCGCACAGGAGCCCGCAGCTTCTCGGCTCGGCCTCGGCC
TCCGCGCGGTGCCCGCAAGGCAGAGTCCCATGTGCAGCCGAGAAGGCGGGCGGGGAGCAGTCC
CCGCTGCTGCACCTCGACCTCTTCAACTTCGACTGCCAGAGGCGGAGGGCAGCCGCTACGTGCTGACC
AGCCCCCGCTCGCTAGAGGCTGCGCCCCGCTGTGCGGTCAAGCCGGTGGAGCTGCTGCCACGGGCCCTG
GCCGACCTGGTGCAGAGGCTCCGGGCCCTCCATGCGGGTGGCCACCGCCTGTATGAGCCTACGAG
GCGGAGCGGCGCAAGCTGCAGCAATGCCGGCCGAGCGGAGCGCATCATGCGCGAGGAGAAGCGG
CGTCTTTTACGCCTTTGAGCCCCGCGGCCGCCGCCGCCGCCGCCGCCGCCGCCGCCGCCGCCGCCGAGC
GCGGGCAGCAGCAGCTGCAGCAGCGCCAGCCTCCCGCCTCGCCCGCCGCGTGGCCCGGCAAG
GCTTCCCTAGTCCCTCCTCGCCCGGACCCAACCTCCGCCAGCGGGTTCTCGAACAGGCAGGAAGAGC
CATTGCTGGACTCACTGTCCCAGCGCGTGGGGCGCCCTCAGCTCCGAGTCGGGCGCATCGTCGTC
TCCTACAGTGGGAAAGCTTGGAGGAGTGCCTGGCCGCTCGGCCAGGAACAGCTGCCCA
GGGGGTGGCGTCTCCACCACCAACGCTCCGGGCCGCCGCTGCCCTGACCTGGTTCCGATACC
GGCCGAGCTTACGCTCGGCGACCTGAGCCATTGCGCCGAGACCGCTCAGCACGTGGAGCGCATCGTG
CGTCAAGTGCCTGAGAGCGCGGCTGCGCGGGGTGCCGGAGCGTACCAGGAAGATCGCGCGCTCATG
CTGGCGCGGCACAGGAGGAGCTCTGCTGCTGGAGCAACGCGCGGCTGCCACAGCCAGTGGGAGCTG
CAGCGCGTGCACGCCAAGCAGCGCGGGAGCGCGAGGAGCGGGAGAAGCAGCGCGCCCTAGAGCAGGGC
CGCCGAGCCTGGCCCGCAGGTGGAGGAGCGCGAGGCCCGCTGGCCGCAAGAGCGCGAGGCGGGC
CGCGCGGAGCGCGGCGAGTACGAGCGCAGCGAGGAGCGCGCGGGAGCTGGCCGAACGCCAGGGCCTC
CTGCGCGGGAGCGGGCGGAGCGCGCGGCCGGGAGGATCGCTGCGCAAGCTTACAGGAGCAGAAC
CTGAAGCAACGTGAGGAAGGTCTACAGGAAGGGCGGGAGCGGGCTGAGCAGATCCGAGGGAGCGCGCC
CAGCGCGCGGCCCGCCAAAGCAGCGGCAGGAGGGCCAGCTGCAGCGGGAGAAGCGGGAGCTGAGCCGG
GCCGAGCGGGCGCCACGAGGCGCTGCTACAGGGTGGACCCGCGCAGCAGCGCCAGGAGCGAGAGGGC
CTGCGGAGCTCGTGGAAGCCAGTTTGGGCCGTGCGCAGGAGAACTACGAGCATTGGTGGAGCAGCGC
ACCCGGGAGTGCGGGAGCGGGCCCGGCGAGAGGAGTGCAGGGTGGCGGGCAAGGAGGCGGAGAG
CGCAAAGAGCGGGAACATCAGGCGCACCTGGAGGCGTGGCCCGGGCGGGGAGCGGCTGCAGCAC
GCGACGCAGTGGCCGAGGAAGCAGTGCAGCAGAAGGCGCGCGCTGGCCAGAGCCGGCTGGAGAAG
GAACGAGCCAGCGGCCAACAAGGAGAAGGTGGAGAGGACGAAGACTGCCCGCGGAGAGCTACTC
CAGGCCATCGGGCGAAGCTGGAGCGCAGCGAGCAGTGCAGCGGAACGGCGCAGTGCCTGGAGAGC
GCCCGCTCCACAGCCCGGCTTCTTCCAGTGGGTGAGAAGGTGCGCGAGGAGACCAACACGCGATCC
TTGACCCGCATGGTGGGGAGGCCAGCTACACGCCAGCTGGACCGCAAATAA
  
```

Restriction Sites: Sgfl-MluI

ACCN: NM_001271507

Insert Size: 2124 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001271507.1](#)

RefSeq Size: 4179 bp

RefSeq ORF: 2124 bp

Locus ID: 56936

UniProt ID: [Q9NQR7](#)

Cytogenetics: 14q24.1

MW: 79.7 kDa