

Product datasheet for **SC333078**

FBXL19 (NM_001271651) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FBXL19 (NM_001271651) Human Untagged Clone
Tag:	Tag Free
Symbol:	FBXL19
Synonyms:	Fbl19; JHDM1C
Vector:	pCMV6-Entry (PS100001)



[View online »](#)

Fully Sequenced ORF: >SC333078 representing NM_001271651.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGGTATGAAAGTCCCCGAAAGGGGGAATCTGGGCCCTCGGCCTTGTGACGCCCCCAATGTCGTCTG
AGCAGCCGGGGGCGGGGGCCGGAGCGCGCCGACGCCGAACCCGCTGCCGCCGCTGCCGGGCCTGTGTG
CGAACTGAGTGCGGGGATTGCCACTTCTGCCGAGACATGAAGAAGTTCGGGGGGCCCCGGGCGCATGAAG
CAGTCGTGCCTGCTCCGGCAGTGCACCTGCCCCCGTGTCTCCACACACAGCTGTGTGCCTCTTGTGTGGG
GAGGCTGGGAAGGAGGACACGGTGGAGGGAGAGGAAGAGAAATTTGGTTTGAGCCTCATGGAGTGTACA
ATCTGCAACGAGATCGTCCACCCCGCTGCCTGAAGATGGGGAAGGCTGAGGGTGTATCAATGCAGAG
ATCCCCAACTGCTGGGAGTGCCTCGCTGCACCCAGGAAGGCCGACAGCAAGGATTAGGTGAGGGG
CCTGGCCGCCGTAGGGCCGACAACGGCGAGGAGGGCGCCAGCTTGGGAGCGGATGGAAGCTGACAGAG
GAGCCACCGCTTCCACCGCCCCCGCCAGGCGCAAGGGCCCCCTGCCTGCCGGGCCCCCCCCGGAGGAC
GTGCCTGGGCCCCCAACGAAAGGAAAGGGAGGAGGGAATGAGCCTCCACCCCAAGGAAAAAGGTG
AAAGGAGGCCGAGAGAGGCACCTGAAGAAGAAACCAAGCCGCTTTGGCTCTGCAGAGGGCCAGCG
GTGCCGTCCCCGTCCCCGAGAGGGAGAAGCTAGAGCGTTTCAAGCGGATGTGCCAGCTGCTGGAACGG
GTGCCTGACACCTCCTCTTCTCCTCGGACTCAGACTCCGACTCCGACTCTTCGGGCACATCGCTGAGT
GAGGACGAAGCCCCCGGCGAGGGCCGGAATGGGCGACGGCCAGCCCGGGGCGAGCTCTGGCGAGAAGGAG
AACCCTGGGGGGCGGGCGGCTGTGCGCCCTGGCAGTGGGGGGCCCTACTCAGCTGGCCCCCTGGGCCCA
GCCCCACCACCCCGGCTCCACAGCTGGAGCGGCACGTGGTGGCGCCCCCGCCTCGAAGCCCTGAGCCC
GACACACTCCCCTTGGCTGCTGGATCCGACACCCCTGCCCCGGGCGCCTGGCTTCGCGTCTTCCAG
CACCTCGGGCCGCGGGAGCTGTGTATCTGCATGCGAGTCTGCCGAATTTGAGCCGCTGGTGCTATGAC
AAGCGTCTGTGGCTCGAATGGACCTGAGCCGGCGGAAGTCACTGACCCCGCCATGCTCAGTGGTGTG
GTTCCGCCAGCCCCGTGCCCTGGACCTCAGCTGGACAGGTGTCTCAAGAAGCAGCTCATGTGGCTT
CTGAACCGACTACAAGGCCTGCAGGAGCTGGTGTCTCTGGCTGCTCCTGGCTCTCTGTCTCTGCCCTG
GGCTCAGCCCCACTGCCAGCCCTGCGGCTCCTGGACCTCCGCTGGATCGAGGATGTTAAAGACTCCAG
CTCCGGGAGTTGCTGCTGCCTCCACCAGACACCAACCAGGGCAAACAGAGAGCCGTGGTGGCTGCAG
GGGGTGGCAGAACTGCTGCTGGCAGTTTGGAGCTGACAGATGCCTCCCTGCGTCTCCTGCTGCGTAC
GCACCCAGCTGAGCGCCCTGGACCTGAGCCACTGCGCCACGTGCGGGACCCAGTGTTACCTCCTC
ACGGCCCCCAGTCCCCACTCCGCGAGACCCTGGTGCACCTCAATCTTGTGGTTGCCACCGCTAACG
GACCACTGCCTCCGCTGTTCCGCCGTGCCCTCGTCTACGCCGCTAGACCTGCGCTCCTGCCGCCAG
CTCTACCCGAAGCTTGTGCCCGCTGGCAGCTGCCGGGCCCCCTGGCCCTTCCGCTGCCCTGAGGAG
AAGCTGCTTCTCAAGGACAGTAG
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001271651

Insert Size: 1956 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001271651.1</u>
RefSeq Size:	3685 bp
RefSeq ORF:	1956 bp
Locus ID:	54620
Cytogenetics:	16p11.2
Protein Families:	Druggable Genome, Transcription Factors
MW:	71.3 kDa
Gene Summary:	This gene encodes a member of the Skp1-Cullin-F-box family of E3 ubiquitin ligases. The encoded protein is reported to bind to the transmembrane receptor interleukin 1 receptor-like 1 and regulate its ubiquitination and degradation. This protein has been linked to the regulation of pulmonary inflammation and psoriasis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2013] Transcript Variant: This variant (2) uses an alternate in-frame splice site compared to variant 1 . The encoded isoform (2) is shorter than isoform 1.