

Product datasheet for **SC333137**

ZNF700 (NM_001271848) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF700 (NM_001271848) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF700
Vector:	pCMV6-Entry (PS100001)



[View online »](#)

Fully Sequenced ORF: >SC333137 representing NM_001271848.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGCCCTGCTGTAGTCACAGGAGCTGTAGAGAGGACCCCGGTACATCTGAAAGCCGGGAAATGATGTTT
CAGGACCCAGTGGCCTTTGAGGATGTGGCTGTGAACCTCACCCAGGAAGAGTGGACATTGCTGGATATT
TCCCAGAAGAATCTCTTCAGGGAAGTGATGCTGGAACCTTCAGGAACCTGACCTCTATAGGAAAAAA
TGGAGTGACCAGAACATTGAATATGAGTACCAAAACCCAGAGAAGCTTCAGGAGTCTCATAGAAGAG
AAAGTCAATGAAATTAAGAAGACAGTCATTGTGGAGAACTTTTACCCAGGTTCCAGATGACAGACTG
AACTTCCAGGAGAAGAAAGCTTCTCCTGAAGTAAATCATGTGACAGCTTTGTGTGTGCAGAAGTTGGC
ATAGGTAACATCATCTTTAATATGAGCATCAGAGGTGACACTGGACACAAGGCATATGAGTATCAGGAA
TATGGACCAAAGCCATATAAGTGTCAACAACCTAAAAATAAGAAAGCCTTCAGGTATCGCCCATCCATT
AGAACACAAGAAAGGGATCACACTGGAGAGAAACCCTATGCTTGTAAAGTCTGTGAAAAACCTTTATT
TTCCATTCAAGCATTGGAAGACACATGGTAATGCACAGTGGGGATGGAAGTTATAATGTAAATTTTGT
GGGAAAGCCTTCCATTCTTTCAGTTTATATCTTATCCATGAAAGAACTCACACTGGAGAGAAACCATAT
GAATGTAAACAATGTGGTAAATCCTTTACTTATTCTGCTACCCTTCAAATACATGAAAGAACTCACACT
GGGAGAAAGCCCTATGAATGTAGCAATGTGATAAAGCATTTCATAGTTCTAGTTCTATCATAGACAT
GAAAGAAGTCACATGGGAGAGAAGCCTTATCAATGCAAAGAATGTGGAAGCATTTCATATACCAAGT
TCTCTTCGTAGACATGAAAGGACCCACTCTGGGAAAAAACCGTATGAATGTAAGCAATATGGGGAAGGC
TTATCCTATCTTATAAGTTTTCAACACACATAAGAATGAACTCTGGAGAAAGACCTTATAAATGTAAG
ATATGTGGGAAAGGCTTTTATTCTGCCAAGTCATTTCAACACATGAAAAAACTCACACTGGAGAGAAA
CGCTATAAATGCAAGCAATGTGGTAAAGCCTTCAATCTTCCAGTTCTTTTCGATATCATGAAAGGATT
CACACTGGAGAGAAACCCCTATGAGTGTAAAGCAGTGTGGGAAAGCCTTCAGATCTGCCTCACAGCTTCA
GTGCACGGTGGGACTCACACTGGAGAGAAACCCTATGAATGTAAGGAATGTGGGAAAGCCTTCAGATCT
ACCTCACACCTTCGAGTGCATGGTAGGACTCATACTGGAGAGAAACCCCTATGAATGTAAGGAATGTGGG
AAAGCCTTCAGATATGTGAAGCACCTTCAAATTCATGAAAGGACAGAAAAACACATAAGAATGCCCTCT
GGAGAAAGACCTTATAAATGTAGTATATGTGAGAAAGGCTTTTATTCTGCCAAGTCATTTCAAACACAT
GAAAAAACTCACACTGGAGAGAAACCCTATGAATGCAACCAATGTGGTAAAGCCTTCAGATGTTGCAAT
TCCCTTCGATATCATGAAAGGACTCACACTGGAGAGAAACCCTATGAGTGTAAAGCAATGTGGGAAAGCC
TTCAGATCTGCCTCACACCTTCGAATGCATGAAAGGACTCACACTGGAGAGAAACCCTATGAGTGTAAAG
CAATGTGGGAAAGCCTTCAGTTGTGCCTCAAACCTTCGAAAGCATGGTAGGACTCACACTGGAGAGAAA
CCCTATGAGTGTAAAGCAATGTGGGAAAGCCTTCAGATCTGCCTCAAACCTTCAGATGCATGAAAGGACT
CACACTGGAGAGAAACCCTATGAATGTAAGGAATGCGAAAAAGCATTCTGTAAATTCTCTTTTCAA
ATACATGAAAGGAAGCACAGAGGAGAGAGGCCCTATGAATGTAAGCATTGTGGGAATGGATTACATCT
GCCAAGATTCTTCAAATACATGCAAGAACACACATTGGAGAGAAACACTATGAATGTAAGGAATGCGGA
AAAGCATTCAATTATTTTCTTCTTGCATATACACGCAAGGACTCATATGGGAGAGAAGCCATATGAA
TGTAAAGATTGTGGGAAAGCATTTCAGTAG
  
```

Restriction Sites: Sgfl-MluI

ACCN: NM_001271848

Insert Size: 2238 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001271848.1</u>
RefSeq Size:	2901 bp
RefSeq ORF:	2238 bp
Locus ID:	90592
UniProt ID:	<u>Q9H0M5</u>
Cytogenetics:	19p13.2
MW:	86.6 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (2) uses an alternate in-frame splice site, compared to variant 1. The encoded isoform (2) is longer than isoform 1.