

Product datasheet for **SC333163**

C18orf10 (TPGS2) (NM_001271949) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: C18orf10 (TPGS2) (NM_001271949) Human Untagged Clone
Tag: Tag Free
Symbol: C18orf10
Synonyms: C18orf10; HMFN0601; L17; PGS2
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC333163 representing NM_001271949.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
ATGGAGGAGGAGGCATCGTCCCCGGGGCTGGGCTGCAGCAAGCCGCACCTGGAGAAGCTGACCCTGGGC
ATCACGCGCATCCTAGAATCTTCCCCAGGTGTGACTGAGGTGACCATCATAGAAAAGCCTCTGCTGAA
CGTCATATGATTTCTTCTGGGAACAAAAGAATAACTGTGTGATGCCTGAAGATGTGAAGAACTTTTAC
CTGATGACCAATGGCTTCCACATGACATGGAGTGTGAAGCTGGATGACCTGGAGGACGATACACATGAA
GCCAGTGTGATCAGCCAGAGAAGCCTCACTTTGACTCTCGCAGTGTGATATTTGAGCTGGATTCATGC
AATGGCAGTGGGAAAGTTTGCCTTGTCTACAAAAGTGGGAAACCAGCATTAGCAGAAGACTGAGATC
TGGTTCTGGACAGAGCCTTATACTGGCATTCTCACAGACACCTTTACTGCCTATTACCGCCTGCTC
ATCACCCACCTGGGCCTGCCAGTGGCAATATGCCTTACCAGCTATGGCATTAGCCCACAGGCCAAG
CAATGGTTCAGCATGTATAAACCTATCACCTACAACACAAACCTGCTCACAGAAGAGACCGACTCCTTT
GTGAATAAGCTAGATCCCAGCAAAGTGTAAAGAGCAAGAACAAGATCGTAATCCCAAAAAGAAAGGG
CCTGTGCAGCCTGCAGGTGGCCAGAAAGGGCCCTCAGGACCCTCCGGTCCCTCCACTTCTCCACTTCT
AAATCCTCCTCTGGCTCTGGAAACCCACCCGGAAGTGA
```

Restriction Sites: Sgfl-Mlul
ACCN: NM_001271949
Insert Size: 798 bp
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001271949.1](#)

RefSeq Size: 3963 bp

RefSeq ORF: 798 bp

Locus ID: 25941

UniProt ID: [Q68CL5](#)

Cytogenetics: 18q12.2

MW: 29.6 kDa

Gene Summary: This gene encodes a protein that is a component of the neuronal polyglutamylase complex, which plays a role in post-translational addition of glutamate residues to C-terminal tubulin tails. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Dec 2012]
Transcript Variant: This variant (2) uses an alternate in-frame splice site in the central coding region compared to variant 1. It encodes isoform 2 which is shorter than isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.