

Product datasheet for **SC333171**

HSP90AB1 (NM_001271970) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	HSP90AB1 (NM_001271970) Human Untagged Clone
Tag:	Tag Free
Symbol:	HSP90AB1
Synonyms:	D6S182; HSP84; HSP90B; HSPC2; HSPCB
Vector:	pCMV6-Entry (PS100001)



[View online »](#)

Fully Sequenced ORF: >SC333171 representing NM_001271970.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGCCTGAGGAAGTGCACCATGGAGAGGAGGAGGTGGAGACTTTTGCCTTTCAGGCAGAAATTGCCCAA
CTCATGTCCCTCATCATCAATACCTTCTATTCCAACAAGGAGATTTTCTTCGGGAGTTGATCTCTAAT
GCTTCTGATGCCTTGGACAAGATTGCTATGAGAGCCTGACAGACCCTTCGAAGTTGGACAGTGGTAAA
GAGCTGAAAATTGACATCATCCCCAACCTCAGGAACGTACCCTGACTTTGGTAGACACAGGCATTGGC
ATGACCAAAGCTGATCTCATAAATAATTTGGGAACCATGCCAAGTCTGGTACTAAAGCATTTCATGGAG
GCTCTTCAGGCTGGTGCAGACATCTCCATGATTGGGCAGTTTGGTGTGGCTTTTATTCTGCCTACTTG
GTGGCAGAGAAAAGTGGTTGTGATCACAAGCACACGATGATGAACAGTATGCTTGGGAGTCTTCTGCT
GGAGGTTCTTCACTGTGCGTGTGACCATGGTGGAGCCATTGGCAGGGGTACCAAAGTATCCTCCAT
CTTAAAGAAGATCAGACAGAGTACCTAGAAGAGAGGGCGGTCAAAGAAGTAGTGAAGAAGCATTCTCAG
TTCATAGGCTATCCCATACCCTTTATTTGGAGAAGGAACGAGAGAAGGAAATAGTGATGATGAGGCA
GAGGAAGAGAAAAGGTGAGAAAGAAGAGGAAGATAAAGATGATGAAGAAAAACCAAGATCGAAGATGTG
GGTTCAGATGAGGAGGATGACAGCGGTAAAGATAAGAAGAAGAAAATAAGAAGATCAAAGAGAAAATAC
ATTGATCAGGAAGAACTAAACAAGACCAAGCCTATTTGGACCAGAAACCTGATGACATCACCACAGAG
GAGTATGGAGAATTCTACAAGAGCCTACTAATGACTGGGAAGACCACTTGGCAGTCAAGCACTTTTCT
GTAGAAGGTCAAGTTGGAATTCAGGGCATTGCTATTTATTCCTCGTCCGGCTCCCTTTGACCTTTTGG
AACAAAGAAGAAAAGAACAACATCAAACCTATGTCCGCCGTGTGTTTCATCATGGACAGCTGTGATGAG
TTGATACCAGAGTATCTCAATTTTATCCGTGGTGTGGTTGACTCTGAGGATCTGCCCTGAACATCTCC
CGAGAAATGCTCCAGCAGAGCAAATCTTGAAGTCATTTCGAAAAACATTGTTAAGAAGTGCCTTGAG
CTCTTCTGAGCTGGCAGAAGACAAGGAGAATTACAAGAAATCTATGAGGCATTCTCTAAAAATCTC
AAGCTTGGAAATCCACGAAGACTCCACTAACCGCCGCCCTGTCTGAGCTGCTGCGTATCATACCTCC
CAGTCTGGAGATGAGATGACATCTGTGACAGTATGTTTCTCGCATGAAGGAGACACAGAAGTCCATC
TATTACATCACTGGTGGAGCAAAGAGCAGGTGGCCAACCTCAGCTTTTGTGGAGCAGTGCAGAAACGG
GGCTTCGAGGTGGTATATATGACCGAGCCATTGACGAGTACTGTGTGCAGCAGCTCAAGGAATTTGAT
GGGAAGAGCCTGGTCTCAGTTACCAAGGAGGGTCTGGAGCTGCCTGAGGATGAGGAGGAGAAGAAGAAG
ATGGAAGAGAGCAAGGCAAAGTTTGAACCTCTGCAAGCTCATGAAAGAAATCTTAGATAAGAAGTT
GAGAAGGTGACAATCTCCAATAGACTTGTGTCTCACCTTGTGCTGATTGTGACCAGCACCTACGGCTGG
ACAGCCAATATGGAGCGGATCATGAAAGCCAGGCACTTCGGGACAACCTCACCATTGGGCTATATGATG
GCCAAAAAGCACCTGGAGATCAACCCTGACCACCCATTGTGGAGACGCTGCGGCAAGGCTGAGGCC
GACAAGAATGATAAGGCAGTTAAGGACCTGGTGGTGTGCTGCTGTTTGAACCGCCCTGCTATCTTCTGGC
TTTTCCCTTGGAGATCCCAGACCCACTCCAACCGCATCTATCGCATGATCAAGCTAGGTCTAGGTATT
GATGAAGATGAAGTGGCAGCAGAGGAACCAATGCTGCAGTTCCTGATGAGATCCCCCTCTCGAGGGC
GATGAGGATGCGTCTCGCATGGAAGAAGTCGATTAG
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001271970

Insert Size: 2175 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001271970.1
RefSeq Size:	2644 bp
RefSeq ORF:	2175 bp
Locus ID:	3326
UniProt ID:	P08238
Cytogenetics:	6p21.1
Protein Families:	Druggable Genome, Stem cell - Pluripotency
Protein Pathways:	Antigen processing and presentation, NOD-like receptor signaling pathway, Pathways in cancer, Progesterone-mediated oocyte maturation, Prostate cancer
MW:	83.3 kDa
Gene Summary:	<p>This gene encodes a member of the heat shock protein 90 family; these proteins are involved in signal transduction, protein folding and degradation and morphological evolution. This gene encodes the constitutive form of the cytosolic 90 kDa heat-shock protein and is thought to play a role in gastric apoptosis and inflammation. Alternative splicing results in multiple transcript variants. Pseudogenes have been identified on multiple chromosomes. [provided by RefSeq, Dec 2012]</p> <p>Transcript Variant: This variant (3) differs in the 5' UTR, compared to variant 1. Variants 1, 2 and 3 encode the same protein (isoform a).</p>