

## Product datasheet for **SC333208**

### EAAT4 (SLC1A6) (NM\_001272088) Human Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** EAAT4 (SLC1A6) (NM\_001272088) Human Untagged Clone  
**Tag:** Tag Free  
**Symbol:** EAAT4  
**Synonyms:** EAAT4  
**Vector:** pCMV6-Entry (PS100001)  
**Fully Sequenced ORF:** >SC333208 representing NM\_001272088.  
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
ATGAGCAGCCATGGCAACAGCCTGTTCCCTGCGGGAGAGCGGCCAGCGGCTGGGCCGGTGGGCTGGCTG
CAGCGGCTGCAGGAAAGCCTGCAGCAGAGAGCACTGCGCACGCGCCTGCGCCTGCAGACCATGACCCCTC
GAGCACGTGCTGCGCTTCTGCGCCGAAACGCCTTCATTCTGCTGACGGTCAGCGCCGTGGTCATTGGG
GTCAGCCTGGCCTTTCCTGCGCCCATATCAGCTCACCTACCGCCAGATCAAGTACTTCTCTTTTCTT
GGAGAGCTTCTGATGAGGATGCTGCAGATGCTGGTGTACCTCTCATTGTCTCCAGCCTGGTCACAGGT
ATGGCATCCCTGGACAACAAGGCCACGGGGCGGATGGGGATGCGGGCAGCTGTGTACTACATGGTGACC
ACCATCATCGCGTCTTCATCGGCATCCTCATGGTCACCATCATCCATCCCGGGAAGGGCTCCAAGGAG
GGGCTGCACCGGGAGGGCCGGATCGAGACCATCCCCACAGCTGATGCCTTCATGGACCTGATCAGAAAT
ATGTTTCCACAAACCTTGTGGAGGCTGCTTCAAACAGTTCAAGACGCAGTACAGCACGAGGGTGGTA
ACCAGGACCATGGTGAGGACAGAGAACGGGTCTGAGCCGGGTGCCTCCATGCCTCCTCCATTCTCAGTG
GAGAACGGAACAGCTTCTGAAAAATGTCACTCGGGCCTTGGGTACCCTGCAGGAGATGCTGAGCTTT
GAGGAGACTGTACCCGTGCCTGGCTCCGCAATGGCATCAACGCCTGGGCTCGTGGTCTTCTCTGTG
GCCTTTGGGCTGGTCATTGGTGGCATGAAACACAAGGGCAGATCCTCAGGGACTTCTCGACAGCCTC
AATGAGGCTATTATGAGGCTGGTGGGCATCATTATCTGGTGA
```

**Restriction Sites:** Sgfl-Mlul  
**ACCN:** NM\_001272088  
**Insert Size:** 939 bp  
**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).  
**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001272088.1](#)

**RefSeq Size:** 1867 bp

**RefSeq ORF:** 939 bp

**Locus ID:** 6511

**UniProt ID:** [P48664](#)

**Cytogenetics:** 19p13.12

**Protein Families:** Transmembrane

**MW:** 34.8 kDa

**Gene Summary:** Sodium-dependent, high-affinity amino acid transporter that mediates the uptake of L-glutamate and also L-aspartate and D-aspartate (PubMed:7791878). Functions as a symporter that transports one amino acid molecule together with two or three Na(+) ions and one proton, in parallel with the counter-transport of one K(+) ion. Mediates Cl(-) flux that is not coupled to amino acid transport; this avoids the accumulation of negative charges due to aspartate and Na(+) symport (By similarity). Plays a redundant role in the rapid removal of released glutamate from the synaptic cleft, which is essential for terminating the postsynaptic action of glutamate (Probable).[UniProtKB/Swiss-Prot Function]  
Transcript Variant: This variant (3) differs in the 5' UTR compared to variant 1. Variants 1 and 3 encode the same isoform (1).