

Product datasheet for **SC333213**

PAQR6 (NM_001272104) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: PAQR6 (NM_001272104) Human Untagged Clone
Tag: Tag Free
Symbol: PAQR6
Synonyms: PRdelta
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC333213 representing NM_001272104.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
ATGCTCAGTCTCAAGCTGCCCAACTTCTCAAGTCCACCAGGTCCCCGGGTGTTCTGGGAAGATGGC
ATCATGTCTGGCTACGCCGCCACCAGCTCGGCTTTGGACTGTGCTCCTCAGCTCCTTCCAGATGACC
AACGAGACGGTCAACATCTGGACTCACTTCTGCCACCTGGTACTTCTGTGGCGGCTCCTGGCGCTG
GCGGGCGGCCCGGCTTCCGTGCGGAGCCGTACCACTGGCCGCTGCTGGTCTTCTGCTGCCCGCTGC
CTTACCCCTTCGCGTCTGCTGCGCGCACACCTTCAGTCCATGTCGCCCCGATGCGCCACATCTGC
TACTTCTCGACTACGGCGGCTCAGCCTCTACAGTCTGGGCTGCGCCTTCCCCTATGCCGCTACTCC
ATGCCGGCTCCTGGCTGCACGGCCACCTGCACCAAGTTCTTTGTGCTGCGCGCCACTCAACTCCTTC
CTGTGCACCGGCTCCTCCTGCTACTCCCGTTTCTGGAGCTGGAAGCCCTGGGCTCAGTAAAGTCTC
CGCACAGGAGCCTTCGCCTATCCATTCTGTTGACAACCTCCCCTCTTTATCGGCTCGGGCTGTGC
TGGGGCAGGGGCCACGGCTGTGGGCAGGAGGCCCTGAGCACCAGCCATGGCTACCATCTCTTCTGCGC
CTGCTCACTGGCTTCTTCTCGCTCCCACCTGCCTGAAAGGCTGGCACCAGGACGCTTTGATTACATC
GGCCACAGCCACCAGTTATCCACATCTGTGCAGTGTGGCACCACCTCCAGCTGGAGGCAGTGTG
GCTGATATGGGATCACGCAGAGCCTGGCTGGCCACACAGGAACCTGCCCTGGGCTGGCAGGCACAGTG
GCCACACTGGTCTTGGCTGCAGCTGGGAACCTACTCATTATTGCTGCTTTCACAGCCACCCTGCTTCGG
GCCCCAGTACATGCCCTCTGCTGCAGGGTGGCCACTGGAGGGGGTACCCAGGCCAAACAACAGTGA
```

Restriction Sites: SgfI-MluI
ACCN: NM_001272104
Insert Size: 1035 bp
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001272104.1</u>
RefSeq Size:	2071 bp
RefSeq ORF:	1035 bp
Locus ID:	79957
UniProt ID:	<u>Q6TCH4</u>
Cytogenetics:	1q22
Protein Families:	Druggable Genome, Transmembrane
MW:	38 kDa
Gene Summary:	<p>Plasma membrane progesterone (P4) receptor coupled to G proteins (PubMed:23763432, PubMed:23161870). Seems to act through a G(s) mediated pathway (PubMed:23161870). Involved in neurosteroid inhibition of apoptosis (PubMed:23161870). May be involved in regulating rapid P4 signaling in the nervous system (PubMed:23763432). Also binds dehydroepiandrosterone (DHEA), pregnanolone, pregnenolone and allopregnanolone (PubMed:23763432, PubMed:23161870).[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (3) uses alternate splice sites at internal exons and initiates translation at an alternate start codon, compared to variant 1. Variants 2 and 3 encode the same isoform (2), which has distinct N- and C-termini and is shorter, compared to isoform 1.</p>