

Product datasheet for **SC333265**

IFTAP (NM_001276722) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: IFTAP (NM_001276722) Human Untagged Clone
Tag: Tag Free
Symbol: IFTAP
Synonyms: C11orf74; HEPIS; NWC
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC333265 representing NM_001276722.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
ATGTCTGCCCATATGTCAGGATTGGAAATAATGGATGAAGATCAATTAATCAAAGACGTCTTGGATAAA
TTCCTTAATTGTCATGAGCAAACATATGATGAAGAATTTCTGAACACTTTTACTCATCTTTCACAAGAG
GATCATGTGTCCAAAAGGGGAGTGTGGAACTGATTCTTCAGAAAACATTTTACCTCAGCAAAGTT
ACTCATAAAAAATGAAGCAGATGACTACCATCTTAGAAAATAAACCATCTTCTTCGTACTTCATCACAA
TGTTTGAAGAACAGGTAGATAATTTCTAGATTTAGAAGATTTGGACATGGATGAAGAGATTAACCC
CAAATGAGTGAGGATTTGCTGCTGCTTCCAGGAGAAGTGGAGCAGGATGTAAGCACCAGCATTCTTCC
TGTATCCCTTTTGTGCCCAGCCTCCTACCTGTGAAGTGAAGCCAAAGCCAGTGTAAAAGAATGGAC
AAACAGACGGAAGAGATACTTGGAGATGAAGTTCAACTTTTTTCACTTGATGAAGAATTTGATTATGAC
AATGTGATGCTAACCTCCAAGTTTAGTCCTGCAGAGATAGAGAACATCAAAGAGCTATGCAAGCAGCAG
AAGAGAAAGGACACCAGCCAGACTTAGAGAAATCCTGTGACTGA
```

Restriction Sites: Sgfl-Mlul
ACCN: NM_001276722
Insert Size: 666 bp
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001276722.1](#)

RefSeq Size: 1080 bp

RefSeq ORF: 666 bp

Locus ID: 119710

UniProt ID: [Q86VG3](#)

Cytogenetics: 11p12

MW: 25.4 kDa

Gene Summary: This gene encodes a protein that was identified as a cellular interacting partner of non-structural protein 10 of the severe acute respiratory syndrome coronavirus (SARS-CoV). The encoded protein may function as a negative regulator of transcription. There is a pseudogene for this gene on chromosome 1. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2013]

Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (a). Variants 1, 2, 3, and 4 encode the same isoform (a). **Sequence Note:** This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.