

Product datasheet for **SC333370**

C6orf48 (NM_001287487) Human Untagged Clone


Product data:

Product Type: Expression Plasmids
 Product Name: C6orf48 (NM_001287487) Human Untagged Clone
 Tag: Tag Free
 Symbol: C6orf48
 Synonyms: D6S57; G8
 Vector: pCMV6-Entry (PS100001)
 Fully Sequenced ORF: >SC333370 representing NM_001287487.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

ATGGAGAGAAGCTTTGTATGGCTGTCATGCTTAGACAGTGATTCCCTGCAACTTGACCTTCAGGCTGGGA
 GAGGTGGAGAGCCATGCCTGTTCTCCTTCTGCTATGGAATTTGCTGACACAATATCTTCCGCCTGGT
 GCTGGGCATATCCTAAGAACTTACAACCTTCTGTATTATCTGTGTGAGCAGCTGCACCTTATTGGG
 GGAAAAATGCCTGAAAATTAG

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001287487

Insert Size: 228 bp



[View online >](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001287487.1</u>
RefSeq Size:	608 bp
RefSeq ORF:	228 bp
Locus ID:	50854
Cytogenetics:	6p21.33
MW:	8.3 kDa