

Product datasheet for SC333505

Protamine 2 (PRM2) (NM_001286359) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Protamine 2 (PRM2) (NM_001286359) Human Untagged Clone
Tag:	Tag Free
Symbol:	Protamine 2
Synonyms:	CT94.2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC333505 representing NM_001286359. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
 ATGGTCCGATACCGCGTGAGGAGCCTGAGCGAACGCTCGCACGAGGTGTACAGGCAGCAGTTGCATGGG
 CAAGAGCAAGGACACCACGGCCAAGAGGAGCAAGGGCTGAGCCCGGAGCACGTCGAGGTCTACGAGAGG
 ACCCATGGCCAGTCTCACTATAGGCGCAGACACTGCTCTCGAAGGAGGCTGCACCGGATCCACAGGCGG
 CAGCATCGCTCCTGCAGAAGGCGCAAAGACGCTCCTGCAGGCACCGGAGGAGGCATCGCAGAGGTCTG
 CCTGCGCCCCCGCCTTGCCTGCATGTCCCTGA
 ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites:	SgfI-MluI
ACCN:	NM_001286359
Insert Size:	309 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).


[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001286359.1](#)

RefSeq Size: 730 bp

RefSeq ORF: 309 bp

Locus ID: 5620

UniProt ID: [P04554](#)

Cytogenetics: 16p13.13

MW: 12.6 kDa

Gene Summary: Protamines substitute for histones in the chromatin of sperm during the haploid phase of spermatogenesis, and are the major DNA-binding proteins in the nucleus of sperm in many vertebrates. They package the sperm DNA into a highly condensed complex in a volume less than 5% of a somatic cell nucleus. Many mammalian species have only one protamine (protamine 1); however, a few species, including human and mouse, have two. This gene encodes protamine 2, which is cleaved to give rise to a family of protamine 2 peptides. Alternatively spliced transcript variants have also been found for this gene. [provided by RefSeq, Sep 2015]
Transcript Variant: This variant (5) uses an alternate donor splice site at the first exon compared to variant 1. The resulting isoform (5) has a distinct C-terminus compared to isoform 1.