

## Product datasheet for **SC333508**

### **C1orf150 (GCSAML) (NM\_001281838) Human Untagged Clone**

#### Product data:

Product Type:	Expression Plasmids
Product Name:	C1orf150 (GCSAML) (NM_001281838) Human Untagged Clone
Tag:	Tag Free
Symbol:	GCSAML
Synonyms:	C1orf150
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC333508 representing NM_001281838. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGACTACATTTGAAAGAAAACTTCAAGATCAAGATAAGAAAAGCCAAGAAGTTTCATCCACTTCTAAT
CAGGAAAACGAGAATGGCAGTGGTTCTGAAGAAGTGTGCTACACTGTCATTAATCACATCCCCATCAG
AGATCCTCCCTGAGCTCCAATGATGATGGCTATGAGAACATTGACTCCCTCACAAGGAAAGTGAGACAG
TTTAGAGAAAGGTCAGAGACAGAATATGCCCTTCTTAGGACTTCTGTTAGTAGGCCTTGTCTCAGCACC
CATGAGCATGATTATGAAGTTGTGTTCCACACTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-Mlul



[View online »](#)

## Plasmid Map:



ACCN: NM\_001281838

Insert Size: 312 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_001281838.1](#)

RefSeq Size: 3918 bp

RefSeq ORF: 312 bp

Locus ID: 148823

Cytogenetics: 1q44

MW: 11.9 kDa

**Gene Summary:** This gene encodes a protein thought to be a signaling molecule associated with germinal centers, the sites of proliferation and differentiation of mature B lymphocytes. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Aug 2013]  
Transcript Variant: This variant (6) contains multiple alternate exons in the 5' UTR, lacks two in-frame exons and initiates translation at a downstream start codon compared to variant 1. The encoded isoform (c) is shorter than isoform a. Variants 4, 5 and 6 encode the same protein. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.