

Product datasheet for **SC333585**

RAG1AP1 (SLC50A1) (NM_001287589) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RAG1AP1 (SLC50A1) (NM_001287589) Human Untagged Clone
Tag:	Tag Free
Symbol:	SLC50A1
Synonyms:	HsSWEET1; RAG1AP1; SCP; slv; SWEET1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC333585 representing NM_001287589. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCGTGGTCTTCACCCTTGGCATGTTCTCCGCCGCCTCTCGGACCTCAGGCACATGCGAATGACCCG
GAGTGTGGACAACGTCCAGTTCTGCCCTTTCTACCACGGAAGTCAAGCTAAGGTGATTCAAATAAA
TCAACCAATGTCTCTCCTACCCACTCACCATTGCTACCCTTCTCACCTCTGCCTCCTGGTGCCTCTAT
GGGTTTCGACTCAGAGATCCCTATATCATGGTGTCCAACTTCCAGGAATCGTCACCAGCTTTATCCGC
TTCTGGCTTTTCTGGAAGTACCCCAAGGAGCAAGACAGGAAGTACTGGCTCCTGCAAACCTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-Mlul



[View online »](#)

Plasmid Map:



ACCN: NM_001287589

Insert Size: 339 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq:	<u>NM_001287589.1</u>
RefSeq Size:	1077 bp
RefSeq ORF:	339 bp
Locus ID:	55974
Cytogenetics:	1q22
Protein Families:	Transmembrane
MW:	13 kDa
Gene Summary:	<p>Mediates sugar transport across membranes. May stimulate V(D)J recombination by the activation of RAG1.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (7) lacks two alternate consecutive exons compared to variant 1, that causes a frameshift. The resulting isoform (g) has a shorter and distinct N-terminus compared to isoform a.</p>