

## Product datasheet for **SC333631**

### C19orf48 (NM\_001290149) Human Untagged Clone


#### Product data:

Product Type: Expression Plasmids  
 Product Name: C19orf48 (NM\_001290149) Human Untagged Clone  
 Tag: Tag Free  
 Symbol: C19orf48  
 Vector: pCMV6-Entry (PS100001)  
 Fully Sequenced ORF: >SC333631 representing NM\_001290149.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
ATGACAGTGTGGAGGCCGTCTTGGAGATCCAGGCCATCACTGGCAGCAGGCTGCTCTCCATGGTGCCA
GGGCCCCGACAGGCCACCAGGCTCATGCTGGGACCCAACCCAGTGCACAAGGACTTGGCTGCTGAGCCAC
ACACCCAGGAGAAGGTGGATAAGTGGGCTACCAAGGGCTTCTGCAGGCTAGGGGAGGAGCCACCCCCG
CTTCCCTATTGTGACCAGGCCTATGGGGAGGAGCTGTCCATACGCCACCGTGAGACCTGGGCCTGGCTC
TCAAGGACAGACCCGCCTGGCCTGGTGTCTCCAGGGTGAAGCAGGCCAGAATCCTGGGGAGCTGCTC
CTGGTTTGA
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM\_001290149

Insert Size: 354 bp



[View online >](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001290149.1</a></u>
<b>RefSeq Size:</b>	2172 bp
<b>RefSeq ORF:</b>	354 bp
<b>Locus ID:</b>	84798
<b>UniProt ID:</b>	<u><a href="#">Q6RUI8</a></u>
<b>Cytogenetics:</b>	19q13.33
<b>MW:</b>	13.1 kDa