

Product datasheet for **SC333717**

C21orf7 (MAP3K7CL) (NM_001286623) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C21orf7 (MAP3K7CL) (NM_001286623) Human Untagged Clone
Tag:	Tag Free
Symbol:	MAP3K7CL
Synonyms:	C21orf7; HC21ORF7; TAK1L; TAKL; TAKL-1; TAKL-2; TAKL-4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC333717 representing NM_001286623. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGTTCAGCTGATTGCACCTTTAGAAGTTATGTGGAACGAGGCAGCAGATCTTAAGCCCCTTGCTCTG
TCACGCAGGCTGGAATGCAGTGGTGAATCATGGCTCACTACAGCCCTGACCTCCTGGGCCAGAGATG
GAGTCTCGCTATTTTCCCAGGTTGGTCTTGAACACCTGGCTTCAAGCAGTCCTCCTGCTTTTGGCTTC
TTGAAGTGCTTGGATTACAGTATTTAGTCTGCAACAAGTTGGCCATGTTGGAGGACAAT
CCAAAGGTCAGCAAGTTGGCTACTGGCGATTGGATGCTCACTCTGAAGCCAAAGTCTATTACTGTGCC
GTGGAATCCCCAGCTCCCTCTGGGTTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-MluI



[View online »](#)

Plasmid Map:



ACCN: NM_001286623

Insert Size: 375 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001286623.1](#)

RefSeq Size: 1582 bp
RefSeq ORF: 375 bp
Locus ID: 56911
Cytogenetics: 21q21.3
Protein Families: Transcription Factors
MW: 13.4 kDa