

Product datasheet for **SC333793**

ZNF302 (NM_001289190) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	ZNF302 (NM_001289190) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF302
Synonyms:	HSD16; MST154; MSTP154; ZNF135L; ZNF140L; ZNF327
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC333793 representing NM_001289190. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTT TAGTGAACCGTCAGAATTTTGT AATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCC GCGATCGCC
ATGTCTCAGGTGACATTTAGTGATGTGGCTATAGACTTCTCTCATGAAGAGTGGGCATGCCTAGATTCT
GCTCAGAGGGACTTATACAAGGATGTGATGGTCCAGAATTATGAGAACCTGGTCTCTGTAGGTCTTTCC
GTAACTAAGCCATATGTGATCATGTTATTGGAGGATGGAAAAGAGCCCTGGATGATGGAGAAAAAACTG
TCAAAAGCTTACCCATTTCTTTATCACACTCTGTTCTGCTTCTGTGAACTTTGGATTCTCTGCTCTA
TTTGAGCATTGTTCAGAAGTCACTGAAATATTTGAGTTGTCAGAACTATGTGTTTCTGGGTGCTTCA
TTCTTATCCAATTCTCCTAATTCCACTGTAGAAGCTTTTTCAAGAAGTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TCAAAGGATGACGACGATAAGGTTAAACGGCCGGC
```

Restriction Sites: Sgfl-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001289190

Insert Size: 396 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001289190.1](#)

RefSeq Size: 2888 bp

RefSeq ORF: 396 bp

Locus ID: 55900

UniProt ID: [Q9NR11](#)

Cytogenetics: 19q13.11

Protein Families: Transcription Factors

MW: 15 kDa

Gene Summary: This gene encodes a member of the zinc-finger protein family. The encoded protein contains seven C2H2-type zinc fingers and a KRAB domain, but its function has yet to be determined. Alternatively spliced transcript variants have been described. [provided by RefSeq, Mar 2014]
Transcript Variant: This variant (12) differs in the 5' UTR, lacks an alternate in-frame exon, and uses an alternate splice junction compared to variant 1. The resulting isoform (e) lacks an alternate internal segment and has a shorter and distinct C-terminus compared to isoform a. Variants 12, 13, 14, and 15 all encode the same isoform (e).