

Product datasheet for **SC333806**

FAM107B (NM_001282703) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FAM107B (NM_001282703) Human Untagged Clone
Tag:	Tag Free
Symbol:	FAM107B
Synonyms:	C10orf45; HITS
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC333806 representing NM_001282703. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCCGAGCCAGACTACATAGAAGATGACAATCCTGAACTCATTAGGCCTCAGAACTGATCAATCCT
GTAAAAACCTCCCGAACCATCAAGATCTTACAGAGAACTTCTTATGAATCAAAAAAGGGTCTTGCT
CCTCAGAACAAACCAGAATTGCAGAAGGTGATGGAAAAAAGAAAACGAGACCAAGTAATAAAGCAGAAG
GAAGAAGAAGCACAGAAGAAGAAATCTGACTTGGAAATAGAGCTATTAACGGCAGCAGAAGTTGGAG
CAGCTTGAAGTGAAGAAGCAGAAATGCAAGAAGCAAGAAAATGCCCCGAGTTTGTGAAGGTGAAA
GGCAATCTCAGGAGAACAGGCCAAGAAGTCGCCAAGCCAGGAGTCTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-MluI



[View online »](#)

Plasmid Map:



ACCN: NM_001282703

Insert Size: 396 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001282703.1](#)

RefSeq Size: 3311 bp
RefSeq ORF: 396 bp
Locus ID: 83641
UniProt ID: [Q9H098](#)
Cytogenetics: 10p13
MW: 15.6 kDa