

## Product datasheet for SC333814

### ZNF302 (NM\_001289189) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF302 (NM_001289189) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF302
Synonyms:	HSD16; MST154; MSTP154; ZNF135L; ZNF140L; ZNF327
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC333814 representing NM_001289189. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTCTCAGGTGACATTTAGTGATGTGGCTATAGACTTCTCTCATGAAGAGTGGGCATGCCTAGATTCT
GCTCAGAGGGACTTATACAAGGATGTGATGGTCCAGAATTATGAGAACCTGGTCTCTGTAGCAGGTCTT
TCCGTAACCTAAGCCATATGTGATCATGTTATTGGAGGATGGAAAAGAGCCCTGGATGAGGAAAAAA
CTGTCAAAGCTTACCCATTTCTTTATCACACTCTGTTCTGCTTCTGTGAACCTTTGGATTCTCTGCT
CTATTTGAGCATTGTTTCAAGTCACTGAAATATTTGAGTTGTCAGAAGTATGTGTTTTCTGGGTGCTT
CATTTCTTATCCAATTCTCCTAATTCCTGCTGTAAGCTTTTTTCAAGAAGTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: Sgfl-MluI



[View online »](#)

## Plasmid Map:



ACCN: NM\_001289189

Insert Size: 399 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_001289189.1](#)

RefSeq Size: 2951 bp

RefSeq ORF: 399 bp

Locus ID: 55900

UniProt ID: [Q9NR11](#)

Cytogenetics: 19q13.11

Protein Families: Transcription Factors

MW: 15.1 kDa

**Gene Summary:** This gene encodes a member of the zinc-finger protein family. The encoded protein contains seven C2H2-type zinc fingers and a KRAB domain, but its function has yet to be determined. Alternatively spliced transcript variants have been described. [provided by RefSeq, Mar 2014]  
Transcript Variant: This variant (11) differs in the 5' UTR, uses an alternate in-frame splice junction, lacks an alternate in-frame exon, and uses another alternate splice junction compared to variant 3. The resulting isoform (d) lacks an alternate internal segment, contains an additional internal aa, and has a shorter and distinct C-terminus compared to isoform a. Variants 10 and 11 both encode the same isoform (d).