

Product datasheet for **SC333844**

UMAD1 (NM_001302348) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	UMAD1 (NM_001302348) Human Untagged Clone
Tag:	Tag Free
Symbol:	UMAD1
Synonyms:	RPA3-AS1; RPA3OS
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC333844 representing NM_001302348. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCC**CGCATCGCC**
 ATGTTTCACTTCTTCAGAAAGCCTCCGGAATCTAAAAAGCCCTCAGTACCAGAGACAGAAGCAGATGGA
 TTCGTCCTTTTAGGAGATACAACAGATGAGCAAAGAATGACAGCAAGAGGCAAACTTCGGACATAGAG
 GCCAACCAACCTTTGGAGACCAACAAAGAAAATTCATCCAGTGTGACTGTATCAGACCCTGAGATGGAA
 AATAAGGCAGGCCAGACTCTGGAGAACAGCTCATTAAATGGCCGAGCTCCTGAGCGATGTGCCGTTCCAC
 CTGGCCCCGCATGTGCTGGCAGTACAGGGCACCATCACTGACCTTCCCGACCACTACTCTCCTATGAT
 GGCAGCGAAAACCTATCACGGTTTTGGTATGATTTCACTCTTGAAAATTCAGTGCTCTGTGATTCA**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001302348
Insert Size:	414 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).


[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001302348.1](#)

RefSeq Size: 2357 bp

RefSeq ORF: 414 bp

Locus ID: 729852

UniProt ID: [C9J7I0](#)

Cytogenetics: 7p21.3

MW: 15.2 kDa