

Product datasheet for **SC333870**

SNX3 (NM_001300929) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SNX3 (NM_001300929) Human Untagged Clone
Tag:	Tag Free
Symbol:	SNX3
Synonyms:	Grd19; MCOPS8; SDP3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC333870 representing NM_001300929. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTGTAAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGCGGAGACCGTGGCTGACACCCGGCGGCTGATCACCAAGCCGCAGAACCTGAATGACGCCTACGGA
CCCCCAGCAACTTCTCGAGATCGATACAAATCTTCTATTTTCAAGCTGAAAGAATCTACTGTTAGA
AGAAGATACAGTGACTTTGAATGGCTGCGAAGTGAATTAGAAAGAGAGCAAGGTCGTAGTTCCCCCG
CTCCCTGGGAAAGCGTTTTTGCCTCAGCTTCTTTTAGAGGAGATGATGGAATATTTGATGACAATTTT
ATTGAGGAAAGAAAACAAGGGCTGGAGCAGTTTATAAACAAGGTCGCTGGTCATCCTCTGGCACAGAAC
GAACGTTGTCTTACATGTTTTACAAGATGAAATAATAGATAAAAGCTATACTCCATCTAAAATAAGA
CATGCCTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-Mlul



[View online »](#)

Plasmid Map:



ACCN: NM_001300929

Insert Size: 423 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001300929.1](#)

RefSeq Size: 1711 bp

RefSeq ORF: 423 bp

Locus ID: 8724

UniProt ID: [O60493](#)

Cytogenetics: 6q21

MW: 16.3 kDa

Gene Summary: This gene encodes a member of the sorting nexin family. Members of this family contain a phox (PX) domain, which is a phosphoinositide binding domain, and are involved in intracellular trafficking. This protein does not contain a coiled coil region, like most family members. This protein interacts with phosphatidylinositol-3-phosphate, and is involved in protein trafficking. A pseudogene of this gene is present on the sex chromosomes. Alternative splicing results in multiple transcript variants encoding distinct isoforms. [provided by RefSeq, Jul 2014]
Transcript Variant: This variant (4) lacks an in-frame portion of the first coding exon, compared to variant 1. It encodes isoform d, which is shorter than isoform a.