

Product datasheet for SC333878

THEG5 (NM_001278575) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	THEG5 (NM_001278575) Human Untagged Clone
Tag:	Tag Free
Symbol:	THEG5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC333878 representing NM_001278575. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
 ATGCTGTACCATACAGTTTTGCCTACAACCTGCTTGTAGCTTCATTGTAGCTGTGAATGCTCCGTTCTCT
 TCCTTTTACCCCTGAAGATTCGCGAACATGCAGGATTGAGCTGTGCGAGGCATAGCATTCCCCACATG
 GTAGACGAAGGCGGATGTGGGTCTGATGGCTCTACTGCCGTCCACCCGGGAGGATGGACGTGTGTCCG
 GGATATCAGAGCTTGGGCTGACCAGCCACATGTCCATGCAGCCCTCCCGCTGGCTGTGTGCACATTT
 TGGGTGTCACTGAAGGATTTTACTGGTGCAGACTTTACAAAGGCAAAATTGTTTCTCAAGGAAGTGGAC
 ATGAAAGATCTAAATCAGGGCTGCTTAAAAAAGCAAGCACAGATTAATGACCTATCTACATCCTACA
 CGGTTTTAG
 ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites:	SgfI-MluI
ACCN:	NM_001278575
Insert Size:	423 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).


[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_001278575.1

RefSeq Size: 1384 bp

RefSeq ORF: 423 bp

Locus ID: 100507527

Cytogenetics: 19q12

MW: 15.3 kDa