

## Product datasheet for **SC333897**

### **C21orf7 (MAP3K7CL) (NM\_001286619) Human Untagged Clone**

#### Product data:


Product Type:	Expression Plasmids
Product Name:	C21orf7 (MAP3K7CL) (NM_001286619) Human Untagged Clone
Tag:	Tag Free
Symbol:	MAP3K7CL
Synonyms:	C21orf7; HC21ORF7; TAK1L; TAKL; TAKL-1; TAKL-2; TAKL-4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC333897 representing NM_001286619. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGATCAGCACAGCCAGGGTACCTGCTGACAAGCCTGTACGCATCGCCTTAGCCTCAATGACGCCTCA
GATGATACACCCCTGAAGACTCCATTCTTTGGTCTTTCCAGAATTAGACCAGCAGCTACAGCCCTG
CCGCCTTGTCATGACTCCGAGGAATCCATGGAGGTGTTCAAACAGCACTGCCAAATAGCAGAAGAATAC
CATGAGGTCAAAAAGGAAATCACCTGCTTGGAGCAAGGAAGAAGGAGCTATTGCCAAGTTAGATCAG
GCAGAAAAGGAGAAGGTGGATGCTGCTGAGCTGGTTCCGGGAATTCGAGGCTCTGACGGAGGAGAATCGG
ACGTTGAGGTTGGCCAGTCTCAATGTGTGGAACAACCTGGAGAACTTCGAATACAGTATCAGAAGAGG
CAGGGCTCGTCTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-Mlul



[View online »](#)

**Plasmid Map:**


**ACCN:** NM\_001286619

**Insert Size:** 429 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001286619.1](#)

RefSeq Size: 2080 bp  
RefSeq ORF: 429 bp  
Locus ID: 56911  
UniProt ID: [P57077](#)  
Cytogenetics: 21q21.3  
Protein Families: Transcription Factors  
MW: 16.4 kDa