

Product datasheet for **SC333915**

CCDC198 (NM_001283060) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CCDC198 (NM_001283060) Human Untagged Clone
Tag:	Tag Free
Symbol:	CCDC198
Synonyms:	C14orf105
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC333915 representing NM_001283060. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GTCGTTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTATTTTGACATCCCCTGGAACACAGAGAAACAAGTATTATTTAAAGGCATCCACCCAAAGACTT
CAAAAGCTTGAACCCATTGACTTGCCACGAGTAATTACTTCAGGAAGACTCCTGAGCCAGCGAGAAGCC
AGGACAATGCACAAAGCAAAGGTAAGAAAAGAAAATGCAAACCTCAATGTATACTTCTGAGAACAGA
CAATATTTGCATAAGATGCAAGTGTGAAATGATCCGTAAAAGACAAGAGGCCCAAATGGAGTTAAAG
AAAAGTCTTCATGGAGAGGCAAGAATTAATAAGCAAAGTCAAGGGACCATAAAGCCAAGAAAACCTT
CAAAGCACCCAAAGGAATGATGACCATGACCTTCTAACCATGTTGCCTGATGAAATCTTGACAGAGGT
CCCGGGCAAACCCACAGTTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-Mlul



[View online »](#)

Plasmid Map:



ACCN: NM_001283060

Insert Size: 435 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001283060.1](#)

RefSeq Size: 1049 bp
RefSeq ORF: 435 bp
Locus ID: 55195
Cytogenetics: 14q22.3
MW: 17 kDa