

Product datasheet for SC333956

ZNHIT3 (NM_001281432) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNHIT3 (NM_001281432) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNHIT3
Synonyms:	Hit1; PEHO; TRIP3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC333956 representing NM_001281432. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGCGTCGCTCAAATGTAGCACCGTCGTGCGTGATCTGCTTGGAGAAGCCAAATACCGCTGTCCA
GCCTGCCGCGTGCCTACTGCTCGGTAGTCTGCTTCCGGAAGCACAAAGAAGCAGTCAACCTGAAACT
CGTCTGTGAGAAAAAATAAGATCAGCTCTTCTACCAAACCGTAAAGCCTGTGAAAAAAGAT
GATGATGACTCTATAGTGATTTTCTCAATAGTGATGAGGAAGAAGACAGAGTTTCTTTGCAGAATTTA
AAGAATTTAGATGGACTCTCGTCTGTACCCAGGCTGGAGTGCAGCGGCGCAATCTCGGCTCACTGCA
ACCTCTCCCTCCCAGATTCAAGCAATTCTCATGCCTCAGCCTCCCGAGCAGCTGGGACTACAGGCGCT
GCCACTACGCCCAACCAATTTTGTATTTTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001281432
Insert Size:	447 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001281432.1](#)

RefSeq Size: 2064 bp

RefSeq ORF: 447 bp

Locus ID: 9326

UniProt ID: [Q15649](#)

Cytogenetics: 17q12

Protein Families: Druggable Genome, Transcription Factors

MW: 16.5 kDa