

Product datasheet for **SC333991**

CCDC182 (NM_001282544) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CCDC182 (NM_001282544) Human Untagged Clone
Tag:	Tag Free
Symbol:	CCDC182
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC333991 representing NM_001282544. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCC**CGCATCGCC**
 ATGGAACCCCTCTACCAGGCTGGGTCCATTCTCATGACGGTGAATACCCTACAGGGGAAAAAATGATA
 GAGAGTGGCCTCCAGTCTGGAGACTTTTCCCTGTCCCAGTCATGGCCCTCCTGCCTCCACCACCCGCT
 GACTTGGAGATCCTGCAGCAGAAGGTGGCCGGGTGCAACGGGAAGTGGAGGACTTTAAGAAAGAGGCA
 TTGAAGTCCATTACCTTGAAGACGCCTTCTGCGAGATGAATGGAGCCCTGGTGCAACAGGAGGAG
 CAGGCGGCTCGGTGAGGCAGCGCTAAGGGAGGAGGAGGACCGTGGCATCGTGCGCAACAAGTTCTC
 ACCTTCTGTGGCGCGGAGAAACAGCTCCGGGAGCACTGCAAGCGGCTGGAGGACCTGCTGCTGGAC
 AGGGGACGTGACGCCCTGCGTGCCACCAAGAAGAGCCAGGCTGAC**TGA**
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites:	SgfI-MluI
ACCN:	NM_001282544
Insert Size:	462 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).


[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_001282544.1

RefSeq Size: 849 bp

RefSeq ORF: 462 bp

Locus ID: 101927581

UniProt ID: A6NF36

Cytogenetics: 17q22

MW: 17.5 kDa