

Product datasheet for **SC334018**

CDV3 (NM_001282763) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CDV3 (NM_001282763) Human Untagged Clone
Tag:	Tag Free
Symbol:	CDV3
Synonyms:	H41
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC334018 representing NM_001282763. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCAAATAAGCAGTGAAAAGGAAGAGACGATAATGAAAAGAGACAAGATCCAGGTGATAACTGGGAA
GAAGGTGGAGGTGGTGGAGGTATGGAAAAATCTTCAGGTCCCTGGAATAAACAGCTCCAGTACAA
GCACCTCCTGCTCCAGTAATTGTTACAGAAACCCAGAACCGATGACTAGTGGTGTGTATAGGCCT
CCTGGGGCCAGGTTAACCACAACAAGGAAAACACCACAAGGACCACCAGAAATCTACAGTGATACACAG
TTCCCATCCCTGCAGTCAACTGCCAAGCATGTAGAAAGCCGGAAGGATAAAGAAATGGAGAAGAGCTTT
GAAGTAGTAAGACACAAAATAGAGGTAGGGATGAGGTTTCAAAAACCCAGGCCCTTAAACTTCAGCTA
GACAACCAATATGCTGTGCTTGAATAATCAGAAAAGCAGCCACTCACAAATCAATTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-Mlul



[View online »](#)

Plasmid Map:



ACCN: NM_001282763

Insert Size: 471 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001282763.1](#)

RefSeq Size: 2899 bp
RefSeq ORF: 471 bp
Locus ID: 55573
UniProt ID: Q9UKY7
Cytogenetics: 3q22.1
MW: 17.3 kDa