

Product datasheet for **SC334025**

CCDC75 (GPATCH11) (NM_001278505) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CCDC75 (GPATCH11) (NM_001278505) Human Untagged Clone
Tag:	Tag Free
Symbol:	CCDC75
Synonyms:	CCDC75; CENP-Y; CENPY
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC334025 representing NM_001278505. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCC**CGGATCGCC**
 ATGTCCAGGAAAAGTGGCATTGGTCATGAGGCATCATTAAACGGAAAGCAGAGGAAAAATTGGAAGC
 TACAGAAAAAGATTACATGAAAAACCAAGCTGAAGAAAAAGCTGCAGAACAGTTTCGAATGCGACTT
 AAAAATAAGCAAGATGAAATGAAGCTAGAAGGAGATCTCAGAAGAAGCCAGCGAGCCTGTCAACAACTG
 GATGTCCAGAAAAATATTCAGGTTCCAGGGAAGCATGGTACTGGTTGAGGCTTGAAGAGGAGACTGAA
 GAAGATGAAGAAGAAAAAGAACAGGATGAAGATGAATATAAGAGTGAAGATTTAAGCGTACTGGAAAAA
 TTACAAATATTGACTAGTTATTTAAGAGAAGAACATCTATATTGTATTTGGTGTGGAACAGCCTATGAA
 GATAAAGAAGACCTATCTTCAAATTGTCCAGGACCAACTCTGCAGATCATGACT**AA**
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites:	SgfI-MluI
ACCN:	NM_001278505
Insert Size:	471 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).


[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_001278505.1

RefSeq Size: 3654 bp

RefSeq ORF: 471 bp

Locus ID: 253635

UniProt ID: Q8N954

Cytogenetics: 2p22.2

MW: 18.6 kDa