

Product datasheet for **SC334026**

KLHL32 (NM_001286254) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	KLHL32 (NM_001286254) Human Untagged Clone
Tag:	Tag Free
Symbol:	KLHL32
Synonyms:	BKLHD5; dj21F7.1; KIAA1900; UG0030H05
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC334026 representing NM_001286254. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTAACCGTCAGAATTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGTGTATGAACCTAACCAAGAATAAGTGGATAAGCCGTAGCCCCATGCTGCAGAGAAGGGTCTACCAT
TCCATGGCTGCTGTACAAAGGAAGCTTTATGTTCTTGGAGGCAATGACCTAGACTACAATAATGACCCG
ATCCTTGTGCGCCATATAGATTCTTACAACATAGACACTGACCAGTGGACACGTTGTAATTTCAACCTG
CTGACTGGCCAAAATGAATCTGGAGTTGCTGTCCATAATGGGAGAATATATTTAGTTGGTGGATTTCA
ATTTGGACAAATGAGCCTCTGGCTTGTATCCAGGTACTGGATGTAAGCAGAGAAGGCAAAGAAGAAGTA
TTCTATGGGCCTACACTCCCTTTTCTTCCAATGGAATAGCAGCATGCTTCCTTCCAGCTCCATATTTT
ACATGCCCTAACCTTCAAACCTTCAAGTGCCTCATCACAGGATTGGCACCATCTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCCGC
```

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001286254
Insert Size:	471 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001286254.2](#)

RefSeq Size: 2765 bp

RefSeq ORF: 471 bp

Locus ID: 114792

UniProt ID: [Q96NJ5](#)

Cytogenetics: 6q16.1

MW: 17.8 kDa