

Product datasheet for **SC334090**

p23 (PTGES3) (NM_001282604) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	p23 (PTGES3) (NM_001282604) Human Untagged Clone
Tag:	Tag Free
Symbol:	PTGES3
Synonyms:	cPGES; P23; TEBP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC334090 representing NM_001282604. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTGTAAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGTCTTTGGAAAAGCAGCCTGCTTCTGCAAAGTGGTACGATCGAAGGGACTATGTCTTCATTGAATTT
TGTGTTGAAGACAGTAAGGATGTTAATGTAATTTTAAAAATCCAAACTTACATTCAGTTGTCTCGGA
GGAAGTGATAATTTTAAGCATTTAAATGAAATTGATCTTTTTCACTGTATTGATCCAAATGATTCCAAG
CATAAAAAGAACGGACAGATCAATTTTATGTTGTTACGAAAAGGAGAATCTGGCCAGTCATGGCCAAGG
TTAACAAAAGAAAGGGCAAAGCTTAATTGGCTTAGTGTGACTTCAATAATTGAAAAGACTGGGAAGAT
GATTCAGATGAAGACATGTCTAATTTTGATCGTTTCTCTGAGATGATGAACAACATGGGTGGTGATGAG
GATGTAGATTTACCAGAAGTAGATGGAGCAGATGATGATTCACAAGACAGTGATGATGAAAAATGCCA
GATCTGGAGTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCC
```

Restriction Sites: Sgfl-Mlul



[View online »](#)

Plasmid Map:



ACCN: NM_001282604

Insert Size: 495 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001282604.1](#)

RefSeq Size: 2094 bp

RefSeq ORF:	495 bp
Locus ID:	10728
UniProt ID:	Q15185
Cytogenetics:	12
Protein Families:	Druggable Genome, Nuclear Hormone Receptor
MW:	19.2 kDa
Gene Summary:	<p>This gene encodes an enzyme that converts prostaglandin endoperoxide H2 (PGH2) to prostaglandin E2 (PGE2). This protein functions as a co-chaperone with heat shock protein 90 (HSP90), localizing to response elements in DNA and disrupting transcriptional activation complexes. Alternative splicing results in multiple transcript variants. There are multiple pseudogenes of this gene on several different chromosomes. [provided by RefSeq, Feb 2016]</p> <p>Transcript Variant: This variant (5) contains an alternate exon in the 5' region, and initiates translation at an alternate start codon, compared to variant 1. The encoded isoform (e) has a distinct and longer N-terminus, compared to isoform a.</p>