

Product datasheet for **SC334097**

CEACAM21 (NM_001288773) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	CEACAM21 (NM_001288773) Human Untagged Clone
Tag:	Tag Free
Symbol:	CEACAM21
Synonyms:	CEACAM3; R29124_1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC334097 representing NM_001288773. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTGTAAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCCCAAGACACTTTCTGTTGGTCACTCTATCCACAGAGTCAGTGGCTCAGCCCTCCATCCAAGCC
AGCAGCACCACAGTCACAGAGAAGGGCTCCGTGGTCTGACCTGCCACACAAATAAACTGGAACCTCT
TTCCAGTGGATTTTCAACAACCAGCGTCTGCAGGTCACGAAGAGGATGAAGCTGTCCTGGTTTAACCAT
GTGCTCACCATAGACCCCATCAGGCAGGAGGACGCTGGGAGTATCAGTGTGAGGTCTCAACCCAGTC
AGCTCAAACAGGAGCGACCCCTCAAGCTGACTGTAAAATCAGATGACAACACTTAGGCATCCTGATC
GGGTCTGGTTGGGAGTCTTCTGGTGGTGCCTTGTGTGTTTCTGCTCCTCCGAAAACTGGCAGG
GCCAGCGATCAGAGTGACTTCAGGGAGCAGCAGCCCCAGCCTCCACCCCGCCATGGACCTCTGAC
AGCTCCATCTCCTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCC
```

Restriction Sites: Sgfl-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001288773

Insert Size: 498 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001288773.2](#)

RefSeq Size: 1711 bp

RefSeq ORF: 498 bp
Locus ID: 90273
UniProt ID: [Q3KPI0](#)
Cytogenetics: 19q13.2
Protein Families: Transmembrane
MW: 18 kDa