

## Product datasheet for SC334121

### ZCCHC17 (NM\_001282572) Human Untagged Clone

#### Product data:


Product Type:	Expression Plasmids
Product Name:	ZCCHC17 (NM_001282572) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZCCHC17
Synonyms:	HSPC251; pNO40; PS1D
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC334121 representing NM_001282572. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAATTCAGGAAGGCCTGAGACCATGGAAAACCTGCCTGCTCTACACTATTTTCCAAGGAGAGGTT
GCTATGGTGACAGACTATGGGGCCTTTATCAAAATCCCAGGCTGTCGGAAGCAAGGTCTGGTCCATCGA
ACTCATATGTCATCCTGTCGGGTGGATAAGCCCTCTGAGATAGTAGATGTTGGAGATAAAGTGTGGGTG
AAGCTTATTGGCCGAGAGATGAAAATGATAGAATAAAAGTATCCCTCTCCATGAAGTTGTCAATCAA
GGGACTGGGAAAGACCTTGATCCCAACAATGTTATCATTGAGAGAAGAAGAAAAAGAAACATAGAGATA
GGAAGTCATCTGACTCTGACAGCTCAGACTCTGAGAGTGATACAGGCAAGAGGGCAAGGCACACATCAA
AAGACAGCAAGGCAGCAAGAAGAAGAAAAAGAAGAAGCACAAGAAGCACAAGGAGTGAGAGT
ATAAAGAGTGTAGGGGTGGTTGA
ACGCGTACGCGGCCGCTCGAGCAGAAAACATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: Sgfl-MluI



[View online »](#)

**Plasmid Map:**


**ACCN:** NM\_001282572

**Insert Size:** 507 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001282572.1](#)

**RefSeq Size:** 1388 bp

RefSeq ORF: 507 bp  
Locus ID: 51538  
UniProt ID: Q9NP64  
Cytogenetics: 1p35.2  
MW: 19.1 kDa