

Product datasheet for **SC334201**

NSL1 (NM_001297739) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	NSL1 (NM_001297739) Human Untagged Clone
Tag:	Tag Free
Symbol:	NSL1
Synonyms:	C1orf48; DC8; MIS14
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001297739, the custom clone sequence may differ by one or more nucleotides

```
ATGGCGGGTCTCCTGAGTTGGTGGTCCTTGACCTCCATGGGACAAGGAGCTCGCGGCTGGCACAGAGA  
GCCAGGCCTTGGTCTCCGCCACTCCCCGAGAAGACTTTTCGGGTGCGCTGCACCTCGAAGCGGGCTGTGAC  
CGAAATGCTACAACGTGCGGCCGCTTCGTGCAAAAGCTCGGGGACGCTCTGCCGAGGAGATTCGGGAG  
CCCCTCTGCGAGATGCGCAGTGGACTTTTGAATCAGCTGTGCAAGAGAATATCAGCATTAAATGGGCAAG  
CATGGCAGGAAGCTTCAGATAATTGTTTTATGGATTCTGACATCAAAGTACTTGAAGATCAGTTTGATGA  
AATCATAGTAGATATAGCCACAAAACGTAAGCAGTATCCAGAAAGATCCTGGAATGTGCATCAAACC  
ATAAAAGCAAAACAAGAAATTCTGAAGCAGTACCACCCTGTTGTACATCCACTGGACCTAAAATATGACC  
CTGATCCAGGAGCAAGCCTGGTGGACTAGCTCTAGGGGCAGTTTAG
```

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001297739
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001297739.1](#), [NP_001284668.1](#)

RefSeq Size: 1573 bp

RefSeq ORF: 537 bp

Locus ID: 25936

Cytogenetics: 1q32.3

Gene Summary: This gene encodes a protein with two coiled-coil domains that localizes to kinetochores, which are chromosome-associated structures that attach to microtubules and mediate chromosome movements during cell division. The encoded protein is part of a conserved protein complex that includes two chromodomain-containing proteins and a component of the outer plate of the kinetochore. This protein complex is proposed to bridge centromeric heterochromatin with the outer kinetochore structure. Multiple transcript variants encoding different isoforms have been found for this gene. There is a pseudogene of the 3' UTR region of this gene on chromosome X. [provided by RefSeq, Jul 2014]

Transcript Variant: This variant (5) lacks a large segment of the 3' coding region and 3' UTR, which results in a frameshift, compared to variant 1. The encoded isoform (4) has a distinct C-terminus and is shorter than isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.