

## Product datasheet for **SC334284**

### ZCCHC10 (NM\_001300817) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ZCCHC10 (NM_001300817) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZCCHC10
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001300817, the custom clone sequence may differ by one or more nucleotides

```
ATGGCGACTCCCATGCATCGGCTAATAGCCCGGAGACAAGCATTGACACAGAGTTGCAGCCTGTCAAGA  
CCTTTTGGATCCTGATTCAGCCATCCATCGTTATTAGTGAAGCAAATAAGCAACATGTAAGATGTCAGAA  
ATGCTTGGAAATTTGGACATTGGACTTATGAATGCACAGGAAAAAGAAAATACCTACATAGGCCCTCAAGG  
ACAGCAGAACTAAAGAAAGCTTTAAAAGAAAAAGAAAACAGATTATTGCAACAAAGCATTGGAGAAA  
CCAATGTAGAAAGAAAGGCCAAGAAAAAAGTTCCAGTAGCAGTAGCAGTGACAGTTCTGCCAGTGATTC  
TTCATCAGAGAGTGAAGAAACATCTACCTCTTCTCCTCAGAGGACAGTGACACTGATGAAAGCTCCTCT  
AGTTCTCATCTCAGCCTCCTCCACAACCTCTTCTCCTCCTCTGATTCAGACTCAGATTCCAGCTCTT  
CCAGTAGCAGTAGCACCAGCACAGATAGCAGCTGACGATGAACCACCAAGAAGAAGAAAAAGAAATA  
G
```

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001300817
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001300817.2](#), [NP\\_001287746.1](#)

**RefSeq Size:** 2240 bp

**RefSeq ORF:** 561 bp

**Locus ID:** 54819

**UniProt ID:** [Q8TBK6](#)

**Cytogenetics:** 5q31.1