

## Product datasheet for **SC334296**

### PRH1 (NM\_001291314) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	PRH1 (NM_001291314) Human Untagged Clone
Tag:	Tag Free
Symbol:	PRH1
Synonyms:	Db-s; PA; PIF-S; Pr1/Pr2; PRH2; PRP-1/PRP-2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001291314, the custom clone sequence may differ by one or more nucleotides

ATGCTTCTGATTCTGCTGTCAGTGGCCCTGCTGGCCTTCAGCTCAGCTCAGGATTTAAATGAAGATGTCA  
 GCCAGGAAGATGTTCCCCTCGTAATATCAGATGGAGGAGACTCTGAGCAGTTCTAGATGAGGAGCGTCA  
 GGGACCACCTTTGGGAGGACAGCAATCTCAACCCTCTGCTGGTGATGGGAACCAGGATGATGGCCCTCAG  
 CAGGGACCACCCAACAAGGAGGCCAGCAGCAACAAGGTCCACCACCTCCTCAGGAAAGCCACAAGGAC  
 CACCCAACAAGGAGGCCAGCAGCAACAAGGTCCACCACCTCCTCAGGAAAGCCACAAGGACCACCCA  
 ACAGGGAGGCCATCCCCCTCCTCTCAAGGAAGGCCACAAGGACCACCCAACAGGGAGGCCATCCCCGT  
 CCTCTCGAGGAAGGCCACAAGGACCACCCAACAGGGAGGCCATCAGCAAGGTCTCCCCACCTCCTC  
 CTGGAAAGCCCCAGGGACCACCTCCCCAAGGGGCCGCCACAAGGACCTCCACAGGGGCAGTCTCCTCA  
 GTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001291314
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).


[View online »](#)

<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u>NM_001291314.1, NP_001278243.1</u>
<b>RefSeq Size:</b>	1174 bp
<b>RefSeq ORF:</b>	564 bp
<b>Locus ID:</b>	5554
<b>UniProt ID:</b>	<u>P02810</u>
<b>Cytogenetics:</b>	12p13.2
<b>Protein Families:</b>	Druggable Genome
<b>Gene Summary:</b>	<p>This gene encodes a member of the heterogeneous family of proline-rich salivary glycoproteins. The encoded preproprotein undergoes proteolytic processing to generate one or more mature isoforms before secretion from the parotid and submandibular/sublingual glands. Multiple distinct alleles of this locus including the parotid isoelectric-focusing variant slow (PIF-s), the parotid acidic protein (Pa), and the double band slow (Db-s) isoforms have been characterized. The reference genome encodes the Db-s allele. Certain alleles of this gene are associated with susceptibility to dental caries. This gene is located in a cluster of closely related salivary proline-rich proteins on chromosome 12. Co-transcription of this gene with adjacent genes has been observed. Alternate splicing of this gene results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Oct 2015]</p> <p>Transcript Variant: This variant (1) represents the longer transcript and encodes the shorter isoform (a, also known as Db-s).</p>