

Product datasheet for **SC334398**

DHRS4 (NM_001282990) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DHRS4 (NM_001282990) Human Untagged Clone
Tag:	Tag Free
Symbol:	DHRS4
Synonyms:	CR; NRDR; PHCR; PSCD; SCAD-SRL; SDR-SRL; SDR25C1; SDR25C2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001282990, the custom clone sequence may differ by one or more nucleotides

```
ATGCACAAGGCGGGGCTGCTAGGCCTCTGTGCCGGGCTTGAATTCGGTGCGGATGGCCAGCTCCGGGA  
TGACCCGCGGGACCCGCTCGCAAATAAGGTGGCCCTGGTAACGGCCTCCACCGACGGGATCGGCTTCGC  
CATCGCCCGCGTTTGGCCAGGACGGGGCCATGTGGTCGTCAGCAGCCGGAAGCAGCAGAATGTGGAC  
CAGGCGGTGGCCACGCTGCAGGGGGAGGGGCTGAGCGTGACGGGCACCGTGTGCCATGTGGGAAGGCGG  
AGGACCGGAGCGGCTGGTGGCCACGACTCTGGACATTAATGTGAAGGCCCCAGCCCTGATGACAAAGGC  
AGTGGTGCCAGAAATGGAGAAACGAGGAGGCGGCTCAGTGGTGATCGTGTCTTCCATAGCAGCCTTCAGT  
CCATCTCTCTCTGGATGGACAAGGAAAAGAGGAAAGCATGAAAGAAACCTGCGGATAAGAAGGTTAG  
GCGAGCCAGAGGATTGTGCTGGCATCGTGTCTTTCCTGTGCTCTGAAGATGCCAGCTACATCACTGGGGA  
AACAGTGGTGGTGGGAGGAACCCCGTCCCGCCTCTGA
```

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001282990
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001282990.1, NP_001269919.1</u>
RefSeq Size:	1078 bp
RefSeq ORF:	600 bp
Locus ID:	10901
UniProt ID:	<u>Q9BTZ2</u>
Cytogenetics:	14q11.2
Protein Families:	Druggable Genome
Protein Pathways:	Metabolic pathways, Retinol metabolism
Gene Summary:	<p>Reduces all-trans-retinal and 9-cis retinal. Can also catalyze the oxidation of all-trans-retinol with NADP as co-factor, but with much lower efficiency. Reduces alkyl phenyl ketones and alpha-dicarbonyl compounds with aromatic rings, such as pyrimidine-4-aldehyde, 3-benzoylpyridine, 4-benzoylpyridine, menadione and 4-hexanoylpyridine. Has no activity towards aliphatic aldehydes and ketones (By similarity).[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (5) lacks two alternate coding exons, but maintains the reading frame, compared to variant 1. The encoded isoform (5) is shorter than isoform 1.</p>