

Product datasheet for **SC334408**

Phospholipase C beta 2 (PLCB2) (NM_001284299) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Phospholipase C beta 2 (PLCB2) (NM_001284299) Human Untagged Clone
Tag:	Tag Free
Symbol:	Phospholipase C beta 2
Synonyms:	PLC-beta-2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001284299, the custom clone sequence may differ by one or more nucleotides

ATGTCTCTGCTCAACCCTGTCCTGCTGCCCCCAAGGTGAAGGCCTATCTGAGCCAAGGGGAGCGCTTCA
 TCAATGGGATGATGAACTACAGTTGCCTCTCCAGTTATCCTCCGTGTGGATCCTAAGGGCTACTACT
 ATACTGGACGTATCAAAGTAAGGAGATGGAGTTTCTGGATATCACCAGCATCCGGGATACTCGTTTGGG
 AAGTTTGCCAAGATGCCAAGAGCCAGAAGCTCCGGGACGTCTTCAACATGGACTTTCCTGATAACAGTT
 TCCTGCTGAAGACACTCACGGTGGTGTCCGGCCCGGACATGGTGGACCTCACCTCCACAACCTTCGTCTC
 CTACAAGGAGAACGTGGCAAGGCCTGGGCTGAGGACGTACTGGCCCTAGTCAAACATCCGCTGACGGCC
 AACGCCTCCCGCAGCACCTTCTGGACAAGATCCTTGTGAAGCTCAAGATGCAGCTCAACTCTGAAGGGA
 AGATTCCGGTGAAGAACTTTTCCAGATGTTTCTGCTGACCGCAAGCGGTGGAAGCTGCTCTCAGTGC
 CTGCCACCTCCCCAAAGGCAAACCTGGAGGAGCGAGATAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001284299
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).


[View online »](#)

Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001284299.1, NP_001271228.1</u>
RefSeq Size:	1840 bp
RefSeq ORF:	600 bp
Locus ID:	5330
Cytogenetics:	15q15.1
Protein Families:	Druggable Genome
Protein Pathways:	Alzheimer's disease, Calcium signaling pathway, Chemokine signaling pathway, Gap junction, GnRH signaling pathway, Huntington's disease, Inositol phosphate metabolism, Long-term depression, Long-term potentiation, Melanogenesis, Metabolic pathways, Phosphatidylinositol signaling system, Taste transduction, Vascular smooth muscle contraction, Wnt signaling pathway
Gene Summary:	<p>The protein encoded by this gene is a phosphodiesterase that catalyzes the hydrolysis of phosphatidylinositol 4,5-bisphosphate to the second messengers inositol 1,4,5-trisphosphate (IP3) and diacylglycerol. The encoded protein is activated by G proteins and has been shown to be involved in the type 2 taste receptor signal transduction pathway. In addition, nuclear factor kappa B can regulate the transcription of this gene, whose protein product is also an important regulator of platelet responses. [provided by RefSeq, Jan 2017]</p> <p>Transcript Variant: This variant (4) lacks multiple 3' coding exons and contains an alternate 3' terminal exon, resulting in a distinct 3' coding region and 3' UTR, compared to variant 1. The encoded isoform (4) is shorter and has a distinct C-terminus, compared to isoform 1.</p>