

## Product datasheet for **SC334450**

### VSTM1 (NM\_001288792) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	VSTM1 (NM_001288792) Human Untagged Clone
Tag:	Tag Free
Symbol:	VSTM1
Synonyms:	SIRL-1; SIRL1; UNQ3033
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001288792, the custom clone sequence may differ by one or more nucleotides

```
ATGACCGCAGAATTCCTCTCCCTGCTTTGCCTCGGGCTGTGTCTGGGCTACGAAGATGAGAAAAAGAATG
AGAAACCGCCCAAGCCCTCCCTCCACGCCTGGCCAGCTCGTGTTGAAGCCGAGAGCAATGTGACCCCT
GAAGTGTGAGGCTCATTCCAGAATGTGACATTTGTGCTGCGCAAGGTGAACGACTCTGGGTACAAGCAG
GAACAGAGCTCGGCAGAAAACGAAGCTGAATCCCTTACGGACCTGAAGCCTAAGGATGCTGGGAGGT
ACTTTTGTGCCTACAAGACAACAGCCTCCCATGAGTGGTCAGAAAGCAGTGAACACTGCAGCTGGTGGT
CACAGATAAACACGATGAACTTGAAGCTCCCTCAATGAAAACAGGTTTCATCATCTGAGGAATCCACCAAG
AGAACCAGCCATTCCAACTTCCGGAGCAGGAGGCTGCCGAGGCAGATTTATCCAATATGAAAAGGGTAT
CTCTCTCGACGGCAGACCCCAAGGAGTGACCTATGCTGAGCTAAGCACCAGCGCCTGTCTGAGGCAGC
TTCAGACACCACCAGGAGCCCCAGGATCTCATGAATATGCGGCACTGAAAGTGTAG
```

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001288792
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001288792.1](#), [NP\\_001275721.1](#)

**RefSeq Size:** 941 bp

**RefSeq ORF:** 618 bp

**Locus ID:** 284415

**UniProt ID:** [Q6UX27](#)

**Cytogenetics:** 19q13.42

**Protein Families:** Transmembrane

**Gene Summary:** Isoform 2: Behaves as a cytokine, promoting IL17A secretion by CD4+ T-cells, and differentiation and activation of IL17 producing helper T-cells (TH17).[UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (3) lacks an in-frame exon in the central coding region, compared to variant 1. The encoded isoform (3) is shorter, compared to isoform 1.