

Product datasheet for **SC334518**

RAG1AP1 (SLC50A1) (NM_001287587) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RAG1AP1 (SLC50A1) (NM_001287587) Human Untagged Clone
Tag:	Tag Free
Symbol:	SLC50A1
Synonyms:	HsSWEET1; RAG1AP1; SCP; slv; SWEET1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001287587, the custom clone sequence may differ by one or more nucleotides

```
ATGCGTGGTCTTCACCCTTGGCATGTTCTCCGCCGGCCTCTCGGACCTCAGGCACATGCGAATGACCCGG
AGTGTGGACAACGTCCAGTTCCTGCCCTTCTCACCACGGAAGTCAAGTGAGGGGCCGACGGCTGCGCAA
CCTGGGCTGGCTGAGTTATGGGGCTTTGAAGGGAGACGGGATCCTCATCGTCAACACAGTGGGTGCT
GCGCTTCAGACCCTGTATATCTTGGCATATCTGCATTACTGCCCTCGGAAGCGTGTGTGCTCCTACAGA
CTGCAACCCTGCTAGGGGCTCTCTCTGGGTTATGGCTACTTTTGGCTCCTGGTACCCAACCCTGAGGC
CCGGCTTCAGCAGTTGGGCTCTTCTGCAGTGTCTTACCATCAGCATGTACCTCTACCACTGGCTGAC
TTGGCTAAGGTGATTCAAATAAATCAACCAATGTCTCTCTACCACTCACCATTGCTACCCTTCTCA
CCTCTGCCTCCTGGTGCCTCTATGGGTTTCGACTCAGAGATCCCTATATCATGGTGTCCAACCTTCCAGG
AATCGTCACCAGCTTTATCCGCTTCTGGCTTTTCTGGAAGTACCCCCAGGAGCAAGACAGGAACTACTGG
CTCCTGCAAACCTGA
```

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001287587
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001287587.1](#), [NP_001274516.1](#)

RefSeq Size: 1383 bp

RefSeq ORF: 645 bp

Locus ID: 55974

UniProt ID: [Q9BRV3](#)

Cytogenetics: 1q22

Protein Families: Transmembrane

Gene Summary: Mediates sugar transport across membranes. May stimulate V(D)J recombination by the activation of RAG1.[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (5) uses an alternate splice junction at the 3' end of the second exon compared to variant 1, that causes a frameshift. The resulting isoform (e) has a shorter and distinct N-terminus compared to isoform a.